



Medical Products Catalogue

Medical consumables, PPE, home care and diagnostic products, sourced from trusted manufacturers and inspected before shipment.

Drainage, Suction & Gastroenterology	Respiratory & Anesthesia
Urology	Hypodermic Injection
Surgery	Gynecological Examination
Wound Dressing	Non-Woven & PPE
Home Care	ECG & Diagnostics

One accountable partner from specification to delivery.

ALIT Medical helps healthcare providers secure reliable medical supplies from trusted manufacturers around the world. With more than a decade of sourcing experience, an established manufacturing network and quality inspection before every shipment, we take responsibility for each order from specification to delivery.

We serve hospitals, government healthcare institutions and healthcare groups across the Americas, Europe and the Middle East. ALIT manages product development, sourcing, factory coordination, quality control, documentation and logistics, so you can focus on patient care.


Global manufacturing partners An established network of manufacturers, selected for each product and market.	Quality inspections Our quality control team inspects every order before it is released for shipment.
Regulatory support Technical and regulatory documentation prepared for your product and destination market.	Made to order Products, sizes and packaging sourced or developed to your requirements.

Product range

01	Drainage, Suction & Gastroenterology	01 / 07
02	Respiratory & Anesthesia	08 / 23
03	Urology	24 / 32
04	Hypodermic Injection	33 / 48
05	Surgery	49 / 56
06	Gynecological Examination	57 / 60
07	Wound Dressing	61 / 74
08	Non-Woven & PPE	75 / 82
09	Home Care	83 / 88
10	ECG	89 / 92
11	Diagnostics	93 / 108

Product codes in this catalogue are ALIT Medical reference codes. Quote them in your enquiry and we will confirm specifications, packaging, lead time and the documentation available for your market.


Drainage Suction Gastroenterology



AL-DR01

Suction Catheter


Connector Type: Cap-cone, Plain type, T-type,
Y-type, Y-type transparent
Size: Fr5 - Fr24, Length: 50cm



AL-DR02

Feeding Tube


Made of 100% medical grade PVC
Size: Fr4 - Fr24



AL-DR03

Nelaton Catheter

Size: Fr6 - Fr24
Length: 40cm for male, 20cm for female



AL-DR04

Rectal Tube

Size: Fr20 - Fr36
Length: 40cm



AL-DR05

Silicone Stomach Tube

Made of 100% medical grade silicone
Size: Fr8 - Fr22
Length: 120cm



AL-DR06

PVC Stomach Tube


Size: Fr6 - Fr26
Length: 120cm



AL-DR07

Duodenal Tube

Closed distal end with four lateral eyes
Size: Fr6 - Fr22



AL-DR08

Umbilical Catheter

Size: Fr4 - Fr8, Length: 50cm
Packed in blister bag



AL-DR15

Suction Connecting Tubing

Size: 3/16", 1/4", 9/32"
Length: 180cm, 300cm, 360cm
and so on



AL-DR16

Suction Yankauer Handle with Connecting Tube

Type: Crown tip, Round tip, Plain type
With vent or without vent
Size: 3/16", 1/4", 9/32"
Length: 180cm, 300cm, 360cm



AL-DR17

Closed Suction Set

Designed for 24 & 72 hours use
Size: Fr6, Fr8, Fr10, Fr12, Fr14, Fr16



AL-DR09

Feeding Tube with Luer Lock Cap

Size: Fr4 - Fr8, Length: 40cm
Size: Fr10- Fr22, Length: 120cm



AL-DR11

Ryle's Tube

With metal ball
Size: Fr6 - Fr22



AL-DR10

Feeding Syringe



AL-DR12

Mucus Extractor

Size: Fr6 - Fr14, Length: 35cm, 40cm, 50cm,
Capacity: 25ml
With / without polo filter membrane



AL-DR13

Mucus Extractor

Size: 40ml, 80ml



AL-DR14

Mucus Extractor

Type: F6, F8, F10, F12,
F14, F16
Capacity: 10ml



AL-DR18

Smooth Suction Yankauer Handle

Type: Crown tip 6mm/8mm,
Plain type 4mm/6mm/8mm
With vent or without vent



AL-DR19

Suction Yankauer Handle

Type: Crown tip, Round tip, Plain type
With vent or without vent



AL-DR20

Wound Drainage System

Size: Available 200ml, 400ml, 800ml (Fr8 - Fr26)
Type: Spring
Packed: Blister bag



AL-DR21

Wound Drainage System

Size: Available 200ml, 400ml (Fr8 - Fr26)
Type: Hollow
Packed: Blister bag



AL-DR22

Silicone Wound Drainage System

Made of 100% medical grade silicone for superior biocompatibility



AL-DR23

Silicone Reservoir

Made of 100% medical grade silicone for superior biocompatibility
Size: 100ml, 150ml, 200ml, 400ml
Sterilized double package



AL-DR27

Chest Drainage Bottle

Size: 1600ml, Sole Chamber



AL-DR28

Chest Drainage Bottle


Size: 2000ml



AL-DR24

Silicone Round Fluted Drainage Tube

Made of 100% medical grade silicone for superior biocompatibility
Size: 10R, 15R, 19R, 24R
Sterilized double package



AL-DR25

Silicone Flat Perforated Drainage Tube

Made of 100% medical grade silicone for superior biocompatibility
Size: 3mm×7mm, 3/4, 80cm,
4mm×10mm, 3/4, 80cm
Sterilized double package



AL-DR29

Chest Drainage System

Disposable high negative pressure suction device is consisted of drainage bottle, vacuum indicator, guiding tube, guiding tube lock plate, concave joint, hanging sling, convex joint, connecting tube, connecting tube lock plate, drainage tube joint, drainage tube, puncture needle And the protective cover is composed.

Description:

1. Continue high negative pressure to fix the edge of the wound and promote tissue growth
2. Reduced postoperative tissue response
3. Reduced the production of hematoma
4. Significantly shorten the exudation period
5. Reduce complications
6. Reduce postoperative fever



AL-DR30

Feeding Bag

Size: 1200ml
Pumping using & Gravity using



AL-DR26

Silicone Round Channel / Fluted Drains

Made of medical grade silicone
Utilize internal channels or fluted, to swiftly remove exudate independent channels facilitate drainage and reduce risk of occlusion made by one time, no connector ensures patients' comfort upon removal radio-opaque line through the length for X-ray visualization
Available with or without "Three Face" stainless steel trocar
Size and length can be customized
Size: Fr10, Fr15, Fr19, Fr24, Fr28
Packing: blister bag



AL-DR31

Suction Liner

Product Characteristic:

Simplicity.
Safety.
Economical.

Indication for Use:

This product is used for collecting the biohazardous fluids which produced in ICU and operation room.

Product Features:


1. Double safety design to ensure no liquid is suctioned back to the tubing and suction unit, effectively prevent cross infection and protect the medical staff safe.
2. It's made of environment- friendly polyethylene. It's safe and non-toxic.
3. Easy to use, can be a single use, also can use with the stander.
4. Available in some specifications: with or without shut off filter, with or without solidifier, 1000ml, 1500ml and 2000ml.

Product Advantages:

1. Soft and transparent (Convenient for observation liquid waste).
2. Double safety design (To prevent fluid into the suction systems by mistake).
3. It can be connected with relevant suction device or vacuum device, soft bag is seal and reliable inside and can be with stand the pressure of 60kg when fully filled with liquid (Prevent infections liquid overflow caused by medical personnel).

Specification

Item	Description
Stand	Four tandem stand 4 canisters (a set)
	Twin tandem stand 2 canisters (a set)
Hanger	Two tandem hanger stand
	Single hanger stand
	Single canister



Respiratory & Anesthesia



AL-RE01

Nasal Oxygen Cannula

Features

Available with different length to suit all patients' requirement available with various types of adult, pediatric, infant and neonate available with a wide selection of prong types the soft curved prong can offer the best possible comfort to patient and the flared type can slow down the flow of oxygen



AL-RE05

Nebulizer with Mouth Piece

Atomization rate is about 0.25/min
Drive gas flow is about 4 to 8 L/min
Atomization particle $\leq 5\mu\text{m}$
More choices of T piece and mouthpiece



AL-RE02

O₂+CO₂ Nasal Cannula

Available with two sizes
Designed for oxygen supply and sampling of exhaled CO₂ gas simultaneously
Soft prongs minimize irritation



AL-RE06

Nebulizer Mask

Atomization rate is about 0.25/min
Drive gas flow is about 4 to 8 L/min
Atomization particle $\leq 5\mu\text{m}$
The standard length of tube is 2.1m, and different length is available
Size: XL, L, M, S



AL-RE03

Oxygen Mask

Features

The turn-up rim can ensure comfortable fit with good seal offered with head strap and adjustable nose clip the standard length of tube is 2.1m, and different length is available


Size: XL, L, M, S



AL-RE07

Adjustable Venturi Mask

Features
Provide various choices to suit the individual patients' requirement
locking ring secures flow setting
concentration for 24%, 26%, 28%, 30%, 35%, 40%, 50%
the standard length of tube is 2.1m, and different length is available
Size: XL, L, M, S



AL-RE04

Non-rebreathing Mask

Features

Provide a high concentration of oxygen offered with head strap and adjustable nose clip the standard length of tube is 2.1m, and different length is available the star lumen tubing can ensure oxygen flow even if the tube is kinked

Size: XL, L, M, S



AL-RE08

Adjustable Venturi Mask

Features
Provide various choices to suit the individual patients' requirements, 24%(blue), 28%(yellow), 31%(white), 35%(green), 40%(pink), 50%(orange), 60%(red) single venturi valve is available with fixed concentration the standard length of tube is 2.1m, and different length is available
Size: XL, L, M, S



AL-RE09

Tracheostomy Mask

Features:
360 degree swivel connector provides maximum patient comfort
Offered with universal connector for suction and oxygen delivery
Size: adult, pediatric



AL-RE10

Oxygen Tubing

Standard, universal and winged
the standard length of tube is 2.1m,
and different length is available



AL-RE11

CPR Mask

It is use for first-aid
For adult or children
Standard connector: 22mm
Hard PP box color: red, blue, green, white, yellow



AL-RE12

CPAP MASK

CPAP is used clinically for the treatment of sleep apnea syndrome(SAS)and related diseases
Three different sizes to choose
Size: S, M, L



AL-RE13

Peak Flow Meter

It is an ideal tool for objective monitoring of asthma patients.



AL-RE15

CONE-SHAPE Connector

Two angle or four angle
Material: ABS



AL-RE14

Humidifier

Effective humidifying performance
provides comfortable breathing
Model: 500ml, 250ml, 100ml



AL-RE16

Nebular Container Mask

It is an auxiliary tool of MDI inhalation treatment,
and it is used to store aerosol
Mask material: silicone or PVC
Size: 160ml adult, pediatric and infant

AL-RE17

Incentive Spirometers

Capacity: 600cc, 900cc, 1200cc, PER/SEC
Product Size: 14×7×14cm
Packing: 50pcs/ctn
Carton Size: 38×38×39.5cm

AL-RE18

Incentive Spirometers

Capacity: 200cc, 400cc, 800cc, 1000cc, 1200cc, PER/SEC
Product Size: 15×9.5×7.5cm
Packing: 50pcs/ctn
Carton Size: 50×31×40cm



AL-RE17

AL-RE18



AL-RE19

PVC Manual Resuscitator

Color: Blue
 Single patient use only
 Built-in 60/40cm H₂O pressure relief valve
 The set includes oxygen reservoir bag, PVC mask and oxygen tubing 100% Medical level PVC material
 Latex-free components
 Additional Accessories (Airway, mouth opener) and customized label are available
 Built-in 30mm exhale port for PEEP, Inlet Valve, or Manometer (Optional)
Size: Adult / Pediatric / Infant
Packing: PP Box / PE Bag / Paper Box



AL-RE20

SEBS Manual Resuscitator

Color: Green
 Single patient use only
 Built-in 60/40cm H₂O pressure relief valve
 The set includes oxygen reservoir bag, PVC mask and oxygen tubing 100% Medical level PVC material
 Latex-free components
 Additional Accessories (Airway, mouth opener) and customized label are available
 Built-in 30mm exhale port for PEEP, Inlet Valve, or Manometer (Optional)
Size: Adult / Pediatric / Infant
Packing: PP Box / PE Bag / Paper Box



AL-RE21

Silicone Manual Resuscitator

Color: Natural
 Autoclave to 134°C helping to prevent cross infection and contamination.
 Built-in 60cm H₂O pressure relief valve
 Latex-free high grade medical silicone material.
 Guaranteed for 5 years or alternatively 20 times autoclave cycles.
 All-in-one intake valve connects to oxygen reservoir bag directly.
 Additional Accessories (Airway, mouth opener) and customized label are available.
 Built-in 30mm exhale port for PEEP, Inlet Valve, or Manometer (Optional)
Size: Adult / Pediatric / Infant
Packing: PP Box / PE Bag / Paper Box



AL-AN01

PVC Anesthesia Mask

Size: 0# white color,
 1# pink color,
 2# yellow color,
 3# green color,
 4# red color,
 5# blue color,
 6# orange color



AL-AN02

Soft Cushion Anesthesia Mask

Size: 0# white color,
 1# pink color,
 2# yellow color,
 3# green color,
 4# red color,
 5# blue color,
 6# orange color



AL-AN03

PVC-free Anesthesia Mask

Size: 1# blue color,
 2# orange color,
 3# pink color,
 4# white color,
 5# yellow color,
 6# green color



AL-AN04

Silicone Anesthesia Mask

Size: 000#, 00#, 0#, 1#, 2#, 3#, 4#, 5#, 6#
 All models are reusable and can be autoclaved repeatedly at 134°C, latex-free



AL-AN05

Catheter Mount

Size: 18cm x Φ15mm / Φ22mm

Type: Normal, Expandable, Reinforced



AL-AN06

Anesthesia Breathing Circuit (Expandable)

Size: Adult, Child

Components: 90° elbow, Screw cap, 3-way connector, Expandable pipe, Tube connector


Tube size: 1.0m / 1.2m / 1.5m / 1.8m



AL-AN11

Continuous positive airway pressure circuits

Length: 1.7m



AL-AN12

Breathing Bag

Size: 0.5L, 1L, 2L, 3L

Latex, Latex-free



AL-AN07

Anesthesia Breathing Circuit (Normal)

Size: Adult, Child

Components: 90° elbow, Screw cap, 3-way connector, Corrugated pipe, Tube connector

Tube size: 1.0m / 1.2m / 1.5m / 1.8m



AL-AN08

Anesthesia Breathing Circuit (Reinforced)

Size: Adult

Components: 90° elbow, Screw cap, 3-way connector, Reinforced pipe, Tube connector

Tube size: 1.5m / 1.8m



AL-AN13

Tracheal HME

Recommended tidal volume(ml): >75

Moisture loss*(mg H₂O/L air): 9.0(@Vt 500ml)

Resistance(pa): 50 at 30 L/min

Bacterial filtration efficiency*(%): /

Viral filtration efficiency*(%): /



AL-AN14

Tracheal HME

Recommended tidal volume(ml): >75

Moisture loss*(mg H₂O/L air): 8.7(@Vt 500ml)

Resistance(pa): 15 at 30 L/min

Bacterial filtration efficiency*(%): /

Viral filtration efficiency*(%): /



AL-AN09

Anesthesia Breathing Circuit with Water trap

Components: 90° elbow, Screw cap, 3-way connector, Stiffle cover, Corrugated pipe, Water trap, Tube connector

Tube size: 1.5m / 1.8m / 2.1m / 2.4m



AL-AN10

Coaxial Breathing Circuit

Size: 1.5m, 1.8m, 2.1m, 2.4m



AL-AN15

Tracheal HME

Recommended tidal volume(ml): >75

Moisture loss*(mg H₂O/L air): 11.3(@Vt 500ml)

Resistance(pa): 10 at 30 L/min

Bacterial filtration efficiency*(%): /

Viral filtration efficiency*(%): /



AL-AN16

Tracheal HME

Recommended tidal volume(ml): >75

Moisture loss*(mg H₂O/L air): 8.9(@Vt 500ml)

Resistance(pa): 20 at 30 L/min

Bacterial filtration efficiency*(%): /

Viral filtration efficiency*(%): /



AL-AN17

Tracheal HME

Recommended tidal volume(ml): >75
Moisture loss*(mg H₂O/L air): 7.8(@Vt 500ml)
Resistance(pa): 20 at 30 L/min
Bacterial filtration efficiency* (%): /
Viral filtration efficiency* (%): /



AL-AN19

HME Filter

Recommended tidal volume(ml): 150-1500
Moisture loss*(mg H₂O/L air): 10.0(@Vt 500ml)
Resistance(pa): 105 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN47

HME Filter

Recommended tidal volume(ml): 150-1500
Moisture loss*(mg H₂O/L air): 9.9(@Vt 500ml)
Resistance(pa): 110 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN48

HME Filter

Recommended tidal volume(ml): 150-1500
Moisture loss*(mg H₂O/L air): 9.2(@Vt 500ml)
Resistance(pa): 100 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN49

Bacterial Filter

Recommended tidal volume(ml): /
Moisture loss*(mg H₂O/L air): /
Resistance(pa): 220 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN18

HME

Recommended tidal volume(ml): 75-1500
Moisture loss*(mg H₂O/L air): 11.2(@Vt 500ml)
Resistance(pa): 30 at 30 L/min
Bacterial filtration efficiency* (%): /
Viral filtration efficiency* (%): /



AL-AN20

HME Filter

Recommended tidal volume(ml): 150-1500
Moisture loss*(mg H₂O/L air): 10.9(@Vt 500ml)
Resistance(pa): 105 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN50

Bacterial Filter

Recommended tidal volume(ml): /
Moisture loss*(mg H₂O/L air): /
Resistance(pa): 70 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN51

Bacterial Filter

Recommended tidal volume(ml): /
Moisture loss*(mg H₂O/L air): /
Resistance(pa): 65 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN52

Bacterial Filter

Recommended tidal volume(ml): /
Moisture loss*(mg H₂O/L air): /
Resistance(pa): 50 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN54

HME Filter

Recommended tidal volume(ml): 75-500
Moisture loss*(mg H₂O/L air): 10.1(@Vt 500ml)
Resistance(pa): 200 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN56

HME Filter

Recommended tidal volume(ml): 15-1500
Moisture loss*(mg H₂O/L air): 10.2(@Vt 500ml)
Resistance(pa): 115 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN58


HME Filter

Recommended tidal volume(ml): 150-1500
Moisture loss*(mg H₂O/L air): 10.0(@Vt 500ml)
Resistance(pa): 130 at 30 L/min
Bacterial filtration efficiency* (%): 99.99
Viral filtration efficiency* (%): 99.99



AL-AN59 Spirometry Filter

Product Specification	Size	
With mouthpiece itself	Machine I.D.	25/28/29.5/30/30.5
	Machine O.D.	30/30.5
Mouthpiece Detachable	Machine I.D.	22/24/25/28/29.5/30/30.5/45.5
	Machine O.D.	26/30/30.5



AL-AN21

Oral/Nasal Endotracheal Tube

With cuff
Size: 2.5mm-10.0mm



AL-AN22

Oral/Nasal Endotracheal Tube

Without cuff
Size: 2.0mm-10.0mm



AL-AN27

Reinforced Endotracheal Tube

Type: With cuff or without cuff
Size: 3.0mm - 9.5mm



AL-AN28

Endotracheal Tube with Pre-loaded Stylet

Uncuffed Size: 2.0mm - 10.0mm
High volume low pressure cuff, Size: 3.0mm - 10.0mm
Low Profile cuff, Size: 2.5mm - 10.0mm



AL-AN23

Nasal Preformed Tracheal Tube

With cuff
Size: 3.0mm-10.0mm



AL-AN24

Nasal Preformed Tracheal Tube

Without cuff
Size: 3.0mm-10.0mm



AL-AN29

Reinforced Endotracheal Tube with Pre-loaded Stylet

Uncuffed, Size: 3.0mm - 10.0mm
High volume, low pressure cuff, Size: 3.0mm - 10.0mm
Low profile cuff, Size: 3.0mm - 10.0mm



AL-AN30

Reinforced Endotracheal Tube with Evacuation Lumen

PU cuff, Size: 6.0mm - 9.0mm
PVC cuff, Size: 6.0mm - 9.0mm



AL-AN25

Oral Preformed Tracheal Tube

With cuff
Size: 3.0mm-10.0mm



AL-AN26

Oral Preformed Tracheal Tube

Without cuff
Size: 2.0mm-10.0mm



AL-AN31


Endotracheal Tube with Evacuation Lumen

Size: 6.0mm-9.0mm




AL-AN32

Endotracheal Tube Holder



AL-AN33

Tracheostomy Tube Holder



AL-AN34

Intubating Stylet

Size: 6# 30.2cm, 10# 38.2cm, 14# 38.2cm



AL-AN35

Endotracheal Tube Introducer (Bougie)

Solid Introducer
Size: 6#: 53.5cm, 10#: 70cm, 15#: 70cm
Hollow Introducer
Size: 10#: 90cm, 15#: 90cm



AL-AN36

Tracheostomy Tube

With cuff
Size: 3.0mm - 10.0mm



AL-AN37

Tracheostomy Tube

Without cuff
Size: 3.0mm-10.0mm



AL-AN38

Oropharyngeal Airway

Type: Guedel
Size: 000#, 40mm, pink color
00#, 50mm, blue color
0#, 60mm, black color
1#, 70mm, white color
2#, 80mm, green color
3#, 90mm, yellow color
4#, 100mm, red color
5#, 110mm, light blue color
6#, 120mm, orange color



AL-AN39

Oropharyngeal Airway

Type: Berman
Size: 40mm, 50mm, 60mm, 70mm,
80mm, 90mm, 100mm, 110mm



AL-AN40

Nasopharyngeal Airway


Made from soft PVC.
Nasopharyngeal airway is designed for easy insertion and maximum patient comfort.
It facilitates suctioning in intensive care, general anaesthesia and also provide an alternative airway during emergency procedures.
Sterilization: EO
Size: I.D3.0MM-I.D9.0MM
Packing: Blister Bag



AL-AN41

Disposable PVC Laryngeal Mask

Size: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#



AL-AN42

Disposable Reinforced PVC Laryngeal Mask

Size: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#



AL-AN43

Disposable Silicone Laryngeal Mask

Size: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#



AL-AN44

Reusable Silicone Laryngeal Mask

Size: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#



AL-AN45

Disposable Reinforced Silicone Laryngeal Mask

Size: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#




AL-AN46

Reusable Reinforced Silicone Laryngeal Mask

Size: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#

Urology






AL-UR01

Urethral Catheter

Made of latex, Silicone coated
One way drainage
Sterilized by EO
For single use only
Available with different length
Size: 6-30FR



AL-UR02

2-way Latex Foley Catheter

Made of latex, Silicone coated
Color coding for size identification, Sterilized by EO
For single use only
Available with Rubber & Plastic valve
Available with 3-5cc, 5-10cc, 5-15cc, 10cc, 15cc, 15-30cc balloon
Available with different types:
Size: 6-30 FR



AL-UR03

3-way Latex Foley Catheter

Made of latex, Silicone coated
Color coding for size identification, Sterilized by EO
For single use only
Available with Rubber & Plastic valve
Available with 3-5cc, 5-10cc, 5-15cc, 10cc, 15cc, 15-30cc balloon
Available with different types:
Size: 16-28FR



AL-UR04

2-way Female Foley Catheter


Made of latex, Silicone coated
Color coding for size identification, Sterilized by EO
For single use only
Available with Rubber & Plastic valve
Available with 3-5cc, 5-10cc, 5-15cc, 10cc, 15cc, 15-30cc balloon
Available with different types:
Size: 12-26FR



AL-UR05

2-way Tiemann Foley Catheter

Made of latex, Silicone coated
Color coding for size identification, Sterilized by EO
For single use only
Available with Rubber & Plastic valve
Available with 3-5cc, 5-10cc, 5-15cc, 10cc, 15cc, 15-30cc balloon
Available with different types:
Size: 12-24FR



AL-UR06

Malecot Catheter

With Funnel end
Made of Siliconized Latex
Wings and/or mushroom shaped tips for catheter retention
Drainage eyes for enhanced passage of fluid or clots
Size: 12-40FR



AL-UR07

Latex T-Drainage Tube


Made of natural latex, silicone coated
Size: Fr10 - Fr26
Packed in individual blister bag



AL-UR08

2-way All Silicone Foley Catheter

100% Made of silicone
Color coding for size identification
Plastic valve only
Sterilized by EO
For single use only
Available with different types:
Size: 6-26 FR



AL-UR09

3-way All Silicone Foley Catheter

100% Made of silicone
Color coding for size identification
Plastic valve only
Sterilized by EO
For single use only
Available with different types:
Size: 16-26FR




AL-UR10

Penrose Drainage Tube

Made of natural latex

Length: 45cm

Size: 1/4", 5/16", 3/8", 1", 7/8", 3/4", 5/8", 1/2"



AL-UR14

Economic Urine Bag

Capacity: 2000ml
without bottom outlet

Adapter with cap

Anti-Reflux Flutter Valve

Inlet tube: O.D 6.5mm; 90cm Length

Sterilization: EO

Packing: 1pc/PE Bag, 10pcs/Inner Bag,
250pcs/Ctn




AL-UR11

All Silicone Male External Catheter

Made of 100% silicone

Size: 25mm, 29mm, 31mm, 36mm, 41mm



AL-UR15

Economic Urine Bag

Capacity: 2000ml

Adapter with cap

Anti-Reflux Flutter Valve

Push Pull Valve

Inlet tube: O.D 6.5mm; 90cm Length

Sterilization: EO

Packing: 1pc/PE Bag, 10pcs/Inner Bag,
250pcs/Ctn



AL-UR12

Male External Catheter

Made of latex

Size: 20mm, 25mm, 30mm, 35mm, 40mm



AL-UR16

Economic Urine Bag

Capacity: 2000ml

Adapter with cap

Anti-Reflux Flutter Valve

T-Tap Valve

Inlet tube: O.D 6.5mm; 90cm Length

Sterilization: EO

Packing: 1pc/PE Bag, 10pcs/Inner Bag,
250pcs/Ctn



AL-UR13

Spigot Catheter(Catheter Plug)

Universal size is molded to fit all sizes of catheters or tubes;

Graduated design made it to easily accommodate

tubing of various internal diameters ;

Ribbed designed provides strong, secure connections;

Ergonomic flanged designed assists withdrawal and maximize grip;

Can be supplied sterile in individual peel pouches,
or even in bulk packaging



AL-UR17

Urinary Drainage Bag

Capacity: 2000ml

Anti-Reflux Flutter Valve

Inlet tube: O.D 9mm; 100cm Length

With ribbon

Twist Turn Valve

Sterilization: EO

Packing: 1pc/PE Bag, 100pcs/Ctn



AL-UR18

Urinary Drainage Bag

Capacity: 2000ml
Inlet tube: O.D 9mm; 100cm Length
 Needleless Sampling Port with Cap,
 with bed sheet clamp
 Anti-Reflux Flutter Valve, Air-Vent on the Bag
 With Single Hanger and Rope Hanger
 Silicone Outlet Tube with Metal Folding Clamp
Sterilization: EO
Packing: 1pc/blister Packing, 40pcs/Ctn



AL-UR22

Urine Meter

500cc Urinometer with Drip Chamber,
 2600ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap,
 with bed sheet clamp, with Slide Clamp
 Universal Double-Hook and Rope Hanger
 T-Tap Outlet Valve
Sterilization: EO
Packing: 1pc/blister Packing, 10pcs/Ctn



AL-UR19

Urinary Drainage Bag

Capacity: 2000ml
Inlet tube: O.D 10mm; 110cm Length
 Needleless Sampling Port with Cap,
 with bed sheet clamp, with Slide Clamp
 Three Parts Anti-Reflux Drip Chamber,
 air-Vent on the Bag
 Universal Double-Hook and Rope Hanger
 T-Tap Outlet Valve, Rigid Dog House
Sterilization: EO
Packing: 1pc/blister Packing, 40pcs/Ctn



AL-UR23

Leg Bag

500ml Leg Bag
 Anti-Reflux Flutter Valve
 Twist Valve
 Two Pre-Assembled Straps
Sterilization: EO
Packing: 1pc/PE Bag, 48pcs/Ctn



AL-UR20

Urinary Drainage Bag

Capacity: 2000ml
Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap,
 with bed sheet clamp, with Slide Clamp
 Anti-Reflux Drip Chamber, Air-Vent on the Bag
 Universal Double-Hook and Rope Hanger
 T-Tap Outlet Valve, Rigid Dog House
Sterilization: EO
Packing: 1pc/blister Packing, 40pcs/Ctn



AL-UR24

Leg Bag

750 ml, non-woven fabric on the back, three cavity
 T valve plug,
Inlet tube: O.D 9mm; 30cm Length
 Needleless sampling port with cap
 Silicone tube 4.5cm
 Two straps are available
Sterilization: EO
Packing: 1pc/blister Packing, 100 pcs/Ctn



AL-UR21

Urine Meter

400cc Two-Chamber Urine meter, 2600ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
 Needleless Sampling Port with Cap,
 with bed sheet clamp, with Slide Clamp
 Universal Double-Hook
 T-Tap Outlet Valve, with plastic cover
Sterilization: EO
Packing: 1pc/blister Packing, 10pcs/Ctn



AL-UR25

Leg Bag

500 ml, 750 ml, 1000 ml
 T valve plug
Inlet tube: O.D 9mm; 30cm Length
 Needleless sampling port with cap
 Two Pre-assembled Straps
Sterilization: EO
Packing: 1pc/blister Packing, 100 pcs/Ctn



AL-UR26

Pediatric Urine Collector

Medical grade PVC material
Size: General Type(Oval shape)100ml / 200ml
 Girl Type(Mickey Mouse Shape)100ml / 200ml
 Anti-Reflux Type 100ml / 200ml
Sterilization: EO
Packing: PE Bag



AL-UR27

Priming Bag

Size: 1000ml, 2000ml

AL-UR28

Colostomy Bag

One-Piece Pouching System

One-Piece Pouching System		Cut-to-Fit 1-1/4" (32mm)	Cut-to-Fit 1-3/4" (44mm)	Cut-to-Fit 2-1/4" (57mm)	Cut-to-Fit 2-3/8" (60mm)	Cut-to-Fit 2-3/4" (70mm)
One-Piece Closed Pouch		—	4002 4003	4001 4004	—	—
		—	4013	—	—	—
One-Piece Drainable Pouch		—	4005	4006	4007	4012
		—	4005V	4006V	4007V	4012V
		—	4005I	4006I	4007I	4012I
		—	4005T	4006T	4007T	4012T
		—	4008	4009	—	—
		—	4008T	4009T	—	—
		4014	4015	—	—	—

Two-Piece Pouching System

Two-Piece Pouching System		Cut-to-Fit 1-1/4" (32mm)	Cut-to-Fit 1-3/4" (44mm)	Cut-to-Fit 2-1/4" (57mm)	Cut-to-Fit 2-3/8" (60mm)	Cut-to-Fit 2-3/4" (70mm)
Two-Piece Closed Pouch		—	4016	4017	—	—
		—	4018	4019	—	—
Two-Piece Drainable Pouch		—	4020	4021	4022	—
		—	4023	4024	4025	—

Two-Piece Barrier

Two-Piece Barrier		Cut-to-Fit 1-1/4" (32mm)	Cut-to-Fit 1-3/4" (44mm)	Cut-to-Fit 2-1/4" (57mm)	Cut-to-Fit 2-3/8" (60mm)	Cut-to-Fit 2-3/4" (70mm)
Hydrocolloid adhesive (Embedded)		—	4026	4027	4028	—

Two-Piece Barrier

Two-Piece Barrier		Cut-to-Fit 1-1/4" (32mm)	Cut-to-Fit 1-3/4" (44mm)	Cut-to-Fit 2-1/4" (57mm)	Cut-to-Fit 2-3/8" (60mm)	Cut-to-Fit 2-3/4" (70mm)
Hydrocolloid Skin Barrier With Non-woven Border (Embedded)		—	4029	4030	4031	—
Hydrocolloid adhesive (Adhesive)		—	4035	4036	4037	—
Convex Skin Barrier		4060	4061	4062	—	—

Urostomy Pouching System

Urostomy Pouching System		Cut-to-Fit 1-1/4" (32mm)	Cut-to-Fit 1-3/4" (44mm)	Cut-to-Fit 2-1/4" (57mm)	Cut-to-Fit 2-3/8" (60mm)	Cut-to-Fit 2-3/4" (70mm)
One-Piece Urostomy Pouch		—	4038	4039	—	—
		—	4040	4041	—	—
Two-Piece Urostomy Pouch		—	4044	4045	—	—

Hypodermic Injection



AL-HY01

Disposable Syringe

Size: 1ml, 3ml, 5ml, 10ml, 20ml, 60ml
With needle or without needle
Luer lock or Luer slip 2 parts or 3 parts



AL-HY02

Tuberculin Syringe

Size: 100ul



AL-HY03

Insulin Syringe

Size: 100ul



AL-HY04

Safety Auto Destruct Syringe

Size: 0.05ml, 0.1ml, 0.5ml



AL-HY05

Beauty Syringe

Material: Medical PC
Luer lock, 1ml



AL-HY06

Auto-disable Syringe with Needle

Size: 0.05 ml, 0.1ml, 0.2ml, 0.3ml, 0.4ml, 0.5ml, 0.6ml, 0.8ml, 1ml, 1.2ml, 2ml, 2.5ml, 3ml, 4ml, 5ml, 6ml, 10ml, 12ml, 20ml



AL-HY07

Safe Clip Syringe

Size: 1ml, 3ml, 5ml, 10ml, 20ml Operation procedure
1. Fix the clip before drawing.
2. Same operation as normal syringe in drawing and injecting
3. After injection, push the clip over the needle and lock it onto the clip base



AL-HY08

Safety Syringe



AL-HY09

Hypodermic Needle


Size: 10-30G



AL-HY10

Irrigation Syringe

Size: 50ml, 60ml



AL-HY15

Disposable Infusion Set

Light Resistance



AL-HY16

Disposable Infusion Set with Regulator

With Y site, With Micro-Regulator
With needle or without needle
With filter or without filter



AL-HY11

Insulin Pen Needle


Size: 30G, 31G, 32G 33G



AL-HY12

Scalp Vein Set

Type: Single & Double wings



AL-HY17

Disposable Blood Transfusion Set

With latex tube
With needle



AL-HY18

Disposable Burette Infusion Set

With burette: 100ml, 150ml
With Y site, With needle



AL-HY13

Disposable Infusion Set

Small chamber, With latex tube
With needle or without needle
With filter or without filter



AL-HY14

Disposable Infusion Set

Small chamber, With Y site
With needle or without needle
With filter or without filter



AL-HY19

Disposable Burette Infusion Set

Light Resistance
With burette: 100ml, 150ml
With Y site, with needle



AL-HY20

Three Way Stopcock

Parts: Body(PC), Plug(PE), Screw cap(ABS)



AL-HY21

Extension Tube

Tube length: 10-150cm



AL-HY22

T-shape Extension Tube

Tube Length: 15cm



AL-HY23

Three Way Stopcock with Extension Tube

Tube Length: 10-150cm



AL-HY28

Combi Stopper

Size: Standard, **Material:** ABS, **Sterile:** EO
Fitted with male luer lock at one for safe and secure connection
and the three way stopcock at the other end
Packing: 1pc/blister, 100pcs/box, 50boxes/ctn

AL-HY29

Screw Cap



AL-HY30

Back Check Valve

Type: male to female, tube to female,
tube to tube, tube to male



AL-HY24

**Extension Line With Positive Pressure
Needle Free Connector**

Type: 1-way, 2-way, 3-way



AL-HY25

Three Way Stopcock

Type: Double platoon, Triple platoon,
Triple platoon with extension tube



AL-HY31 I.V. Start Kit

1pc IV cannula dressing
1pc TPE tourniquet
1roll PE tape
1pc Gauze swabs
1pc Povidone-iodine prep pad



AL-HY32

Tourniquet

Type: Plastic cuff
Size: Adult, Child



AL-HY26

Needle Free Connector

Type: Straight, Y type, Intubation



AL-HY27

Heparin Cap

Type: A: yellow, B: Transparent, C: Octagonal
Material: PC



AL-HY33 Latex-free Tourniquet

Thickness: 0.06mm
Width: 3/4", 1"
Length: 18"
Color: Blue, yellow, green



AL-HY34

Pre-cut Roll Tourniquet

Material: TPE



AL-HY35

Spinal Needle

Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
Length: 90mm



AL-HY36

I.V. Cannula

Type: with injection port and suitable wings
Size: 14G-26G, PTEF / FEP / PU flex catheter,
Sterilization: EO
Packing: Hard Blister Packing



AL-HY37

I.V. Cannula

Type: without injection port and with suitable wings.
Size: 14G-26G, PTEF / FEP / PU flex catheter,
Sterilization: EO
Packing: Hard Blister Packing



AL-HY38

I.V. Cannula

Type: Pen type
Size: 14G-26G, PTEF / FEP / PU flex catheter,
Sterilization: EO
Packing: Hard Blister Packing



AL-HY39

Latex Tube

Size: I.D.: 3.0-8.0, O.D.: 5.0-12.0
Length: 20m, 25m



AL-HY40

Epidural Needle

#12-#18(18G-15G)L65-150mm
Needle material: Stainless steel
Needle hub: transparent PP

Color coded hub for size identification Packed in individual Blister Bag
1000pcs/ctn, **Dimension:** 42×27×52cm

AL-HY41

Epidural Catheter

AL-HY42

Connect Adapter Luer Lock



AL-HY43

Anesthetic Products Kit

Standard Kits Include:

1. Drainage Catheter
2. Guide-wire
3. Vessel Dilator
4. Introducer Needle
5. Introducer Syringe
6. Injection Needle
7. Injection Cap
8. Connecting tube

Optional Compound Kits Include:

1. Drainage Catheter
- Standard Kit Accessories
2. 5ml Syringe
3. Surgical Gloves
4. Surgical Pledget
5. Surgical Sheet
6. Surgical Towel
7. Sterile Brush
8. Gauze Pad
9. Suture of Needle
10. Wound Dressing
11. Scalpel



AL-HY44

Blood Purification Products Kit

Standard Kits Include:

1. Hemodialysis Catheter
2. Vessel Dilator
3. Introducer Needle
4. Syringe
5. Guide-wire
6. Adhesive Wound Dressings
7. Heparin Caps
8. Scalpel
9. Needle with Suture

Optional Compound Kits include:

1. Hemodialysis Catheter
- Standard Kit Accessories
2. 5ml Syringe
3. Surgical Gloves
4. Surgical Pledget
5. Surgical Sheet
6. Surgical Towel
7. Sterile Brush
8. Gauze Pad



AL-HY45

Blood Purification Products Kit

Characteristics:

1. The Catheter is made from Premium silicone, which has excellent biocompatibility.
2. The dacron cuff can protect immigration of catheter and reduce the chance of CRBI(Catheter Related Blood stream Infections).
3. The silicone material makes catheter more floppy, hence it can minimize the trauma to the vein and reduce the incident of clot.

Standard Kits Include:

1. Silicone Catheter
2. Vessel Dilator
3. Valved Peelable Introducer
4. Introducer Needle
5. Guide-wire
6. Tunneler (Tocar)
7. Adhesive Wound Dressings
8. Heparin caps
9. Gauze Sponges
10. Scalpel
11. Syringe
12. Needle with Suture



Vacuum Blood Collection Tube

REMARKS:

1. THE INTERIOR OF THE TUBE IS COATED WITH SILICONE.
2. CAN DO ACCORDING TO YOUR SPECIAL REQUEST AND PACKING.



ALP14-1001

Name: Plain Tube
Cap Color: Red

Description:

1. There is no any additive (no anticoagulation) in the tube, which is also called no additive tube. Red closure.
2. **Specimen:** serum. The tube is used for clinical biochemistry and immunology.
3. **Specs:** 13x75 mm; 13x100 mm; 16x100 mm
4. **Draw volume:** 3-10 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±10%
7. Centrifuge for 5-10 minutes after blood coagulation.



ALP14-1002

Name: Glucose Tube
Cap Color: Grey

Description:

1. The tube is using Sodium Fluoride and Potassium Oxalate as additive. Grey closure.
2. **Specimen:** plasma. It's used for glucose determinations.
3. **Specs:** 13x75 mm
4. **Draw volume:** 2 ml, 3 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±10%
7. Invert the tube 6-8 times gently.



ALP14-1003

Name: Procoagulation Tube
Cap Color: Orange / Red

Description:

1. The tube is coated with micronized particles as clot activator to accelerate clotting. Orange closure.
2. **Specimen:** serum. It is used for serum determinations in chemistry, serology, and immunohematology.
3. **Specs:** 13x75 mm; 13x100 mm; 16x100 mm
4. **Draw volume:** 3-10 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±10%
7. Invert the tube 6-8 times immediately after drawing to mix the blood and clot activator thoroughly.
8. Centrifuge for 5-10 minutes after 10-30 minutes' coagulation.



ALP14-1004

Name: ESR Tube
Cap Color: Black

Description:

1. The tube is using 3.8% (0.129 mol/L) Tri-Sodium Citrate as anticoagulation. Black closure.
2. **Specimen:** whole blood. It's used for Westergren ESR (Erythrocyte Sedimentation Rate) test.
3. **Specs:** 13x75 mm; 8x120 mm (glass)
4. **Draw volume:** 1.6, 2.4, 4 ml; 1.28, 1.6 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±4%
7. The mix ratio is 1 part citrate to 4 parts blood, and must draw required volume.
8. Invert the tube 6-8 times gently.



ALP14-1005

Name: PT Tube
Cap Color: Blue

Description:

1. The tube is using 3.2% (0.109 mol/L) Tri-Sodium Citrate as anticoagulation. Blue closure.
2. **Specimen:** whole blood or plasma. It's used for routine coagulation testing.
3. **Specs:** 13x75 mm
4. **Draw volume:** 1.8, 2.7, 3.6 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±4%
7. The mix ratio is 1 part citrate to 9 parts blood, and must draw required volume.
8. Invert the tube 6-8 times gently.



ALP14-1006

Name: Heparin Tube
Cap Color: Green

Description:

1. The tube is using Sodium Heparin or Lithium Heparin as additive. Green closure.
2. **Specimen:** plasma. It's used for plasma determinations in chemistry.
3. **Specs:** 13x75 mm; 13x100 mm; 16x100 mm
4. **Draw volume:** 3-10 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±10%
7. Draw required volume and invert the tube 6-8 times gently.



ALP14-1007

Name: EDTA Tube
Cap Color: Purple

Description:

1. The tube is using EDTAK2 or EDTAK3 as additive. Lavender closure.
2. **Specimen:** whole blood. It's used for hematology determinations.
3. **Specs:** 10x60 mm (glass); 13x75 mm; 16x100 mm
4. **Draw volume:** 1 ml; 2-10 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±10%
7. Invert the tube 6-8 times immediately after drawing.



ALP14-1008

Name: GEL&CLOT Activator Tube
Cap Color: Yellow

Description:

1. The interior of the tube is coated with clot activator (silica particles) to accelerate clotting. A barrier polymer-gel is present at the tube bottom, which will form a barrier separating serum from fibrin and cells during centrifugation. It can prevent substance exchange between cell and serum, guaranteeing biochemical character and chemical components of serum unchanged for a long time.
2. **Specimen:** serum. It's used for serum determinations in chemistry.
3. **Specs:** 13x75 mm; 13x100 mm; 16x100 mm
4. **Draw volume:** 3-9 ml
5. **Storage temperature:** 4-25°C
6. **Max. vacuum degree tolerance:** ±10%
7. Invert the tube for 5-8 times and centrifuge for 5-10 minutes at 1100-1500g after 10-30 minutes' coagulation.



ALP14-1009



STP14-1009A



ALP14-1010



STP14-1010A



STP14-1011



STP14-1012

CAT NO.	NAME	DESCRIPTION
ALP14-1009	Multi-Sample Needle	The uniquely sharpened faceting simplifies penetration of tissue. In comparison with other common needles, painful trauma is considerably reduced. The specially sharpened point results in minimal tissue damage and thromboplastin seepage.
ALP14-1009A	Flash-back Blood Needle	Intuitive, easy-to-use, single-handed activation device combines safety with added reassurance of flashback for the confidence of success.
ALP14-1010	Blood Collection Sets	Blood Collection Sets allows you to deliver exceptional patient care without compromising personal safety. This product is intended for blood collection and infusion procedures. It is a necessity for automatic blood collection in the future.
ALP14-1010A	Venous Blood Specimen Collection Needle	Blood Collection Sets allows you to deliver exceptional patient care without compromising personal safety. This product is intended for blood collection and infusion procedures. It is a necessity for automatic blood collection in the future.
ALP14-1011	Needle Holder	Needle Holder is used with collection needle. It helps to reduce the pain and ensure the safety when collecting the blood sample from patient. As the needle is fixed with needle holder, it will reduce the potential danger when the needle is used on the patient.
ALP14-1012	Safety Needle Holder	The Safety Needle Holder is an improved product of needle holder. It's safe and easy to use. It helps to reduce the pain and ensure the safety when collecting the blood sample from patient.



Flat Bottom



CAT NO.	TUBE NAME	CAP COLOR	TUBE SIZE (MM)	MATERIAL	ADDITIVE	DRAW VOLUME
ALP14-1011PP	PLAIN TUBE	WHITE	15X45	PP	NO ADDITIVE	2.5ML
ALP14-1012PP	EDTA TUBE	RED	15X45	PP	EDTA K2, EDTA K3	2.5ML
ALP14-1013PP	GLUCOSE TUBE	YELLOW	15X45	PP	POTASSIUM OXALATE+ SODIUM FLUORIDE	2.5ML
ALP14-1014PP	HEPARIN TUBE	BLUE	15X45	PP	SODIUM HEPARIN/ LITHIUM HEPARIN	2.5ML



AL-DE01

Dental Needle


Type: Straight, Curved
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 27G
Packing: 1pc/blister bag, 100pcs/box



AL-DE02

Dental Needle

Size: 27G, 30G
Packing: 1pc/blister bag, 100pcs/box



AL-DE07

Salive Ejector

Material: PVC,
Length: 15cm
Packing: 100pcs/bag



AL-DE08

Dental Syringe

Material: PP
Curved Tip: 12ml



AL-DE03

Disposable Bib

Material: 2×20g paper laminated 20g PE film, etc.
Size: 33×46cm, 33×66cm, 40×45cm, etc.
Color: white, green, pink, blue, purple, etc.
Packing: 125pcs/bag



AL-DE04

Dental Tray

Packing: 1set/plastic bag, sterile with E.O gas



AL-DE09

Disposable Dental Mirror (Glass)

Packing Details: (single)100pcs/box, 20boxes/ctn
(bulk)100pcs/bag, 20bags/ctn
CTN Size: (single)50×39×24cm
(bulk)42×30×28cm
GW.: (single)14.2kgs
(bulk)13.5kgs



AL-DE14

Surgical Aspirator Tips

Color: blue, yellow, white
Size: 148×Φ11mm
Material: plastic
Packing Details: 100pcs/bag, 20bags/ctn
CTN Size: 42×33×31cm
GW.: 8kgs



AL-DE05

Dental Sponges

Cotton filled, **Mesh:** 19×9, 20×12, etc.
Size: 5×5cm-8ply, 7.5×7.5cm-8ply, 10×10cm-8ply, etc.
Sterile or non-sterile
Packing: 100pcs/pack



AL-DE06

Dental Cotton Rolls


Material: 100% cotton, natural, super absorbent
A: 20g/bundle, **Dia:** 7cm,
B: 25g/bundle, **Dia:** 7.2cm, etc



AL-DE18

Denture Box

Color: purple, blue, green
Size: 84×75×78mm
Material: PP
Packing Details: 180cases/ctn, 144cases/ctn
CTN Size: 63×40×43cm, 56×36×53cm
GW.: 14kgs, 12kgs



AL-DE21

Dental Impression Tray

Color: Colors
Material: PC/ABS
Packing Details: 10/bag, 100bags/ctn
CTN Size: 63×41×44cm
GW.: 15.1kgs

Material: PP

Description: This product is a scleroid plastic container which is used for collecting the needles of the syringes, infusion sets and pointed instruments like blades, scalp vein sets, suture needle, small glass products, etc. This product is for single use and the medical waste put in it is not allowed to be taken out and must be recycled by the waste disposal unit of the hospital, and then safely and completely burned within 48 hours. The product has the characteristics of non-toxic, safe, cast-resistant, and non-leaky.

Sharp Container



AL-SBX01

Capacity: 4L
Carton Size: 53x42x45cm
Packing: 24pcs/ctn



AL-SBX02

Capacity: 3L
Carton Size: 47x32x45cm
Packing: 30pcs/ctn




AL-SBX03

Capacity: 1.5L
Carton Size: 38x33x48cm
Packing: 60pcs/ctn



AL-SBX04

Capacity: 0.8L
Carton Size: 38x38x39.5cm
Packing: 90pcs/ctn



AL-SBX05

Capacity: 3L
Carton Size: 36x36x37cm
Packing: 30pcs/ctn



AL-SBX06

Capacity: 2L
Carton Size: 45x30x34cm
Packing: 36pcs/ctn



AL-SBX07

Capacity: 1L
Carton Size: 47x36x57cm
Packing: 100pcs/ctn



AL-SBX08

Capacity: 0.7L
Carton Size: 56x24x43cm
Packing: 100pcs/ctn



AL-SBX09

Capacity: 2L
Carton Size: 58x29x29cm
Packing: 40pcs/ctn



AL-SBX10

Capacity: 3L
Carton Size: 47x50x40cm
Packing: 54pcs/ctn



AL-SBX11

Capacity: 4.6L
Carton Size: 58x56x28cm
Packing: 20pcs/ctn



AL-SBX12

Capacity: 2L
Carton Size: 41x34x34cm
Packing: 30pcs/ctn



AL-SBX13

Capacity: 3L
Carton Size: 41x34x31.5cm
Packing: 30pcs/ctn



AL-SBX14

Capacity: 6L
Carton Size: 51x42x34.5cm
Packing: 24pcs/ctn



AL-SBX15

Capacity: 5L
Carton Size: 51x42x25.5cm
Packing: 24pcs/ctn




AL-SBX16

Capacity: 23L
Carton Size: 48x39x65cm
Packing: 20pcs/ctn



AL-SBX17A

Capacity: 10L
Carton Size: 52x25x73cm
Packing: 12pcs/ctn



AL-SBX17B

Capacity: 10L
Carton Size: 51x42x54cm
Packing: 24pcs/ctn



AL-SBX18

Capacity: 8L
Carton Size: 44x38x44cm
Packing: 24pcs/ctn



AL-SBX19

Paper Safety Box
Capacity: 5L
Carton Size: 62x30x258.5cm
Packing: 75pcs/ctn
G.W./N.W.: 20/18.5kgs



AL-SBX20

Needle Destroyer without blade 11x9.5cm
Needle Destroyer with blade 11x9.5cm



AL-SBX21

Vomit Bag
Material: PP
Size: 1500ml



AL-SBX22

Specimen Bag

Color: Red / Yellow / Black
Material: HDPE / LDPE, virgin or recycled
Thickness: 40-150 microns per side
Packaging: Roll packed, or loosed packed
Primary Packaging: labels, polybags
Secondary Packaging: Corrugated cardboard outer cartons, or PP woven bags
Type: Open with adhesive tape, 61x61cm, 61x76cm, 61x80cm, 96.5x120cm
Plain top without adhesive tape, 42x48+9cm, 19x36+5.5cm

AL-SBX23

Specimen Bag

Material: VIRGIN LDPE
Regular Size: 5.5"x6.5", 6"x9", 10"x10", or customized sizes
Thickness: 40-100micron per side
Advantages: 1. Superior quality and competitive price, waterproof and recyclable material.
2. It is environmental protection and excellent quality
3. very tough, Good stretchability, Flexibility,
4. Timely delivery.
5. Each inquiry will be consulted carefully.
6. Easy to open and take.
Packaging: By normal corrugated cardboard box, about the details as customize.
Usually 100pcs/bag, 10bags/ctn or 20bags/ctn



Surgery



AL-SU01

Disposable Surgical Blade

Size: 10#, 10A#, 11#, 12#, 15#, 20#, 21#, 22#, 23#, 24#
Material: Carbon steel or stainless steel
Packing: 1pc/aluminum foil bag, 100pcs or 200pcs/box



AL-SU02

Skin Graft Blade

Size: 158mm
Material: Carbon steel or stainless steel
Packing: 1pc/aluminum foil bag, 100pcs /box



AL-SU03

Stitch Cutter Blade

Size: 65mm, 110mm
Material: Carbon steel or stainless steel
Packing: 1pc/aluminum foil bag, 100pcs/box



AL-SU04

Disposable Scalpel

Size: 10#, 36#
Material: Carbon steel or stainless steel
Packing: 1pc/aluminum foil bag, 10pcs/box



AL-SU05

Safety Surgical Scalpel

Size: 10#, 36#
Material: Carbon steel or stainless steel
Packing: 1pc/aluminum foil bag, 10pcs/box



AL-SU06

Surgical Handle

Size: 3# 12.1cm, 4# 14cm
Material: Stainless steel
Packing: 1pc/aluminum foil bag



AL-SU07

Safety Lancet

Ergonomic design offers a more secure hold when to activate the device
Tri-bevel needle tip ensures the sharpness and accuracy
The spring buffer the puncture making the procedure a more comfortable experience
The complete range of needle sizes from 18G to 32G ensure the accurate
Blood samples for all popular meters
The needle retracts immediately after use
Sterilization: Gamma
Size: 21G, 23G, 26G



AL-SU12

Twist Lancets

Conventional twist-off cap design
Sterilized by gamma-radiation
Smooth tri-bevel needle tip for comfortable sampling experience
Versatile with most lancing devices
Size: 21G, 23G, 26G, 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn

AL-SU13

Flat Lancets (Exclusive to Roche Diagnostics)

5 twist lancets (28G) are used with Roche Accu-Chek Softclix Lancing Device and the needle tip is siliconized for easy and smooth motion during lancing
Size: 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn



AL-SU14

Disposable Electrosurgical Pencil

170mm-ABS pencil appearance
7 cm stainless steel tip
2.36mm-diameter of the tip
3m-PVC medical cable
Maximal output current: 1A
Frequency confine: 0-1.0MHZ
1pc/blister bag



AL-SU15

Surgical Hand Brush

With iodine or without iodine
Packing: 1/mold bag,
20 bags/box, 10boxes/ctn



AL-SU16

Surgical Hand Brush

Material: PP+Nylon



AL-SU08

Pressure Activated Safety Lancets

Needle is fully shielded before and after use
Pre-loaded, pressure activated design
Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization: Gamma
Size: 21G, 23G, 26G



AL-SU09

Button Activated Safety Lancets

Needle tip is sealed with cap through injection moulding
Pre-loaded, button operated design
Single-use only, sterilised by gamma irradiation, cannot be reused
High quality, ultra sharp tri-bevel needle, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization: Gamma
Size: 21G, 23G, 26G



AL-SU10

Lancing Device (Releasable)

AL-SU11

Lancing Device

The user can install the lancet with the device easily. And the device will work with the lancet smoothly without any blockage.
The lancet won't fall on its own when to fire the device, Retracting automatically after firing, the launcher is safe and reliable.



Thread Length: 45cm, 60cm, 75cm, 90cm, 150cm etc.
Thread Diameter: 8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2, 3 etc.
Needle Length: 6mm, 8mm, 12mm, 18mm, 35mm, 40mm, 50mm etc.
Needle Curvature: Straight, 1/2circle, 1/2circle(double), 1/4circle, 1/4circle(double), 3/8circle, 3/8circle(double), 5/8circle, loop round
Cross Section: round bodied, round bodied(heavy), curved cutting, curved cutting(heavy), reverse cutting, reverse cutting(heavy) etc.

Surgical Suture

- AL-SU17 Chromic Catgut
- AL-SU18 Plain Catgut
- AL-SU19 Polyglycolic Acid
- AL-SU20 Silk braided
- AL-SU21 Nylon monofilament
- AL-SU22 Polyester braided
- AL-SU23 Polypropylene



NEEDLE SHAPES

	POINT TYPE	SYMBOL
	taper point	●
	blunt taper point	●
	cutting edge	▲
	reverse cutting edge	▼
	tapercut	⊗
	micro-point spatula curved	◩

1/2 CIRCLE




3/8 CIRCLE



1/4 CIRCLE



5/8 CIRCLE



STRAIGHT



AL-SU17
Chromic Catgut



AL-SU18
Plain Catgut



AL-SU19
Polyglycolic Acid



AL-SU20
Silk braided



AL-SU21
Nylon monofilament



AL-SU22
Polyester braided



AL-SU23
Polypropylene



AL-SU24
Disposable Medical Razor

Material: Carbon steel or stainless steel
Type: Single or double blade



AL-SU25
Disposable Medical Razor

Material: Carbon steel or stainless steel
Type: Single or double blade



AL-SU26

Surgical Glove

Powder or powder-free

Size: 6mm, 6.5mm, 7mm, 7.5mm, 8mm, 8.5mm, 9mm

Packing: 50pairs/box, 500pairs/ctn



AL-SU27

Latex Examination Glove

Powder or powder-free

Size: S, M, L

Packing: 100pcs/box, 1000pcs/ctn



AL-SU28

Vinyl Glove

Size: S, M, L, XL

Packing: 100pcs/box, 1000pcs/ctn



AL-SU29

Gynecology Examination Glove

Material: Natural latex

Size: S, M, L

Packing: paper-paper bag



AL-SU30

Nitrile Examination Glove

Size: S, M, L, XL

Packing: 100pcs/box, 1000pcs/ctn



AL-SU31

TPE Examination Glove

Size: S, M, L, XL

Packing: 100pcs/box, 1000pcs/ctn



AL-SU32

EVA Examination Glove

1.2gm/pc, white color

Package: 1pc/bag, 100bag/inner bag



AL-SU33

PE Glove

Material: Recycled or virgin polyethylene

Thickness: 0.018mm, 0.02mm, 0.022mm

Size: S, M, L

Type: HDPE or LDPE, plain or embossed



AL-SU34

PE Two Fingers Cots

1.2gm/pc, white color

Package: 1pc/bag, 100bag/inner bag

Gynecological Examination

Vaginal Speculum



AL-GE01 Side Screw Type

Size: S, M, L



AL-GE02 Fastener Type

Size: S, M, L



AL-GE03 Spanish Type

Size: S, M, L



AL-GE04 Middle Screw Type

Size: S, M, L



AL-GE05 Push Pull Type

Size: S, M, L



AL-GE06 American Type

Size: S, M, L



AL-GE07

Cervical Brush




AL-GE08

Cervical Brush



AL-GE09

Cervical Brush



AL-GE10

Cervical Brush



AL-GE11

Cervical Brush



AL-GE12

Cervical Scraper




AL-GE13

Cervical Scraper



AL-GE14

Cervical Scraper



AL-GE15

Wooden Tongue Depressor

Size: Adult (150mm×18mm×1.6mm)
Child (140mm×14mm×1.6mm)
EO Sterile / Non-Sterile
Packing: Single medical paper wrapped/Bulk
100pcs/Box, 5000pcs/Ctn



AL-GE16

Cervical Scraper (Wooden)

EO Sterile / Non-Sterile
Packing: Single medical paper wrapped/Bulk
100pcs/Box, 5000pcs/Ctn



AL-GE17

Gynecological Kit

Configuration:
1pc speculum
1pc Gynaecology pad(40×60cm)
1pair latex glove
1pc cervical spatula



AL-GE18

Condom

Specification:
Surface: smooth, dotted, ribbed lubricated
Size: L, M, S (55mm±2, 52mm±2, 49mm±2).
Be available with spermicide nonoxynol 9.
Be available with various fruit smell, such as apple/ banana/ strawberry/ cherry etc.
Be available with square or rectangle packing :
If quantity is large enough, the condom is available with various color such as natural/ red/pink/green etc.
Specification is available with En600, ISO4074, ASTM, WHO, SABS, etc.



AL-GE19

Breast Reliever

Material: PVC+rubber



AL-GE20

Umbilical Cord Clamp Cutter

Material: Pom + Stainless Steel Blade
Feature: Disposable Medical Use
Colour: White, Green, Blue
Function: To cut off the umbilical cord Clamp when a woman giving birth to a baby



AL-GE21

Umbilical Cord Clamp Cutter

Material: ABS
Length: 14cm



AL-GE22

Umbilical Cord Clamp

Size: 5.0 / 5.5 / 6.0cm
Material: ABS/PP/PE
Color: Blue, White, Pink



AL-GE23

Amniotomy Hook

Material: ABS
Type: Straight, Curved.



AL-IVD07

HCG Test Strip

Size: 2.5mm, 3.0mm, etc.
Response Time: 5 minutes
Sensitivity: 25mIU/ml
Specificity: 99%
Packing: 1pc/aluminum foil bag



AL-IVD08

HCG Test Device

Size: 3.0mm, 3.5mm, etc.
Response Time: 5 minutes
Sensitivity: 25mIU/ml
Specificity: 99%
Packing: 1pc/aluminium foil bag



AL-IVD09

HCG Midstream

Size: 4.0mm, 4.5mm, 5.0mm, etc.
Response Time: 5 minutes
Sensitivity: 25mIU/ml
Specificity: 99%
Packing: 1pc/aluminium foil bag


Wound Dressing



AL-WD01

Gauze Swabs

Mesh: 19×15, 26×18, 30×20, etc.
Size: 5cm×5cm, 7.5cm×7.5cm, 10×10cm, etc.
Layer: 8ply, 12ply, 24ply, etc.
Folded or non folded edge
Packing: Sterile or non sterile



AL-WD02

Non-woven Swabs

Size: 5cm×5cm, 7.5cm×7.5cm, 10×10cm, etc.
Layer: 4ply, 6ply, etc.
Packing: Sterile or non sterile



AL-WD03

Trach Sponges

Material: 100% cotton
Mesh: 19×11, 19×15, 20×12, 24×20, etc.
Size: 5cm×5cm, 7.5cm×7.5cm, 10cm×10cm, etc.
Layer: 4ply, 8ply, 12ply, 16ply, etc.
Edge: Folded or unfolded edge
Packing: Sterile or non-sterile



AL-WD04

Cutting Gauze

Mesh: 19×15, 26×18, 30×20, etc.
Size: 10cm×10cm, 15cm×15cm, 20cm×20cm, etc.
Packing: Non sterile



AL-WD05

Gauze Roll

Mesh: 19×15, 26×18, 30×20, etc.
Size: 90cm×91m or 90cm×100m, etc.
Packing: individual blue kraft paper



AL-WD06

Gauze Bandage

Mesh: 19×15, 26×18, 30×20, etc.
Size: 5cm×6y, 7.5cm×6y, 10cm×6y, etc.
Packing: white paper per roll, blue kraft paper per dozen



AL-WD07

Lap Sponge

Mesh: 26×18, 24×20, etc.
Size: 45×45cm-4ply, 30×30cm-4ply, 10×45cm-4ply, etc.
Washed or non-washed, With Blue loop or without, with X-ray or without
Packing: Sterile or non-sterile



AL-WD08

Cotton Roll

Size: 50g, 100g, 200g, 454g, 500g, 1000g, etc.
Packing: Individual blue kraft paper



AL-WD09

Gauze Ball

Mesh: 30×20, 26×18, etc.
Size: 10cm×10cm, 15cm×15cm, 20cm×20cm, 30cm×30cm, etc.
With or without x-ray



AL-WD11

Zig Zag Cotton

Size: 50g, 100g, 200g, 454g, 500g, 1000g, etc.
Packing: Individual Plastic bag



AL-WD10

Cotton Ball

Weight: 0.3g, 0.5g, 1.0g, 2.0g, 3.0g, 4.0g, 5.0g, etc.
Sterile or non-sterile



AL-WD12

Cotton Tip Applicators

Material of the stick: PP, PS, Wooden
Single tip for medical use



AL-WD13

Triangular Bandage

Material: 100% cotton or non-woven
Mesh: 40×34, 50×30, 48×48, etc.
Size: 36"×36"×51", 40"×40"×56", etc.
Packing: 1pc/sealed poly bag



AL-WD14

Thermal Accident Blanket

Size: 132cm×210cm, 140cm×220cm, 160cm×210cm, etc.
Color: Gold, Silver



AL-WD15

O.R. Towel

Constructed of absorbent cotton fiber
Color: Blue, white, green
All the products have different types: Pre-sterile, sterile by steam and by Gamma and non-sterile



AL-WD16

Gauze Eye Pad

Type: Self-adhesive or without adhesive
Size: 5cm×6cm, 8cm×6cm, etc.
Packing: 1pc/pouch, 50pouch/inner box



AL-WD21

Cotton Crepe Bandage

Material: 100% Cotton, **Color:** natural, flesh or bleach
Weight: 75g/m², 80g/m², 85g/m², 90g/m², etc.
Size: 5cm×4m, 7.5cm×4m, 10cm×4m, etc.
Packing: 1 roll / poly bag



AL-WD22

High Elastic Bandage

Material: Polyester + Rubber string
Color: natural, bleach, **Weight:** 90g/m², 105g/m²
Size: 5cm×4m, 7.5cm×4m, 10cm×4m, etc.
Packing: 1 roll / poly bag



AL-WD17

Non-woven Eye Pad

Type: Self-adhesive or without adhesive
Size: 5cm×6cm, 8cm×6cm, etc.
Packing: 1pc/pouch, 50pouch/inner box



AL-WD18

Cosmetic Cotton Swabs

Shape: Round or Square
Type: Plain, Lattice and flower pattern, etc.
Size: Ø5.5cm, Ø6.0cm, 6cm×6cm, etc.



AL-WD23 Net Bandage

Material: Polyester+Rubber string, **Color:** white
Weight: 120g/m², **Size:** 0.5#: 0.8cm×25m, 1#: 1.7cm×25m, 2#: 2cm×25m, 3#: 2.2cm×25m, 4#: 2.6cm×25m, 5#: 3cm×25m, 5.5#: 3.8cm×25m, 5.8#: 4.4cm×25m, 6#: 5.5cm×25m, 7#: 6.8cm×25m, 8#: 8cm×25m, etc. **Packing:** 1 roll / poly bag



AL-WD24

Tubular Bandage

Material: 100% cotton, **Color:** Natural white
Size: 4.5cm×10m, 6.25cm×10m, 6.5cm×10m, etc.
Packing: 1 roll / poly bag



AL-WD19

Elastic Crepe Bandage

Material: Cotton + Spandex
Color: natural, flesh or bleach
Weight: 75g/m², 80g/m², 85g/m², 90g/m², etc.
Size: 5cm×4m, 7.5cm×4m, 10cm×4m, etc.
Packing: 1 roll / poly bag



AL-WD20


PBT Elastic Bandage

Material: Polyester
Color: Bleach
Weight: 28g/m², 30g/m², 60g/m², etc.
Size: 5cm×4m, 7.5cm×4m, 10cm×4m, etc.
Packing: 1 roll / poly bag



AL-WD25 Self-adhesive Elastic Bandage

Material: Cotton with Spandex, Crepe or plain
Color: white, skin
Size: 5cm×4m, 7.5cm×4m, 10cm×4m, etc.
Packing: 1 roll / poly bag



AL-WD26

First-aid Dressing

Size: 6cm×8cm, 8cm×10cm, 10cm×12cm, etc.
Packing: Sterile package or vacuum sterile package



AL-WD27

Non-woven Self-adhesive Bandage

Material: Non-woven PP
Color: white, blue, green, red, brown, yellow, etc.
Size: 5cmx4m, 7.5cmx4m, 10cmx4m, etc.
Packing: 1 roll / poly bag



AL-WD31

Strip Sports Adhesive Tape

Size: 5cmx5m, 7.5cmx5m, 10cmx5m, etc.



AL-WD28

Cotton Self-adhesive Bandage

Material: Cotton
Color: white, blue, green, red, brown, yellow, etc.
Size: 5cmx4m, 7.5cmx4m, 10cmx4m, etc.
Packing: 1 roll / poly bag



AL-WD32

Muscle Tape

Material: Elastic fabric
Color: white, skin color, black, red, green, etc.
Size: 2.5cmx5m, 5cmx5m, 7.5cmx5m, 10cmx5m, etc.
Packing: in hanging box



AL-WD29

Plaster of Paris Bandage

Material: Polyester
Color: Bleach
Weight: 28g/m², 30g/m², 60g/m², etc.
Size: 5cmx4m, 7.5cmx4m, 10cmx4m, etc.
Packing: 1 roll / poly bag



AL-WD33

Zinc Oxide Adhesive Plaster

Material: 10% zinc oxide with 90% hot melt adhesive
Color: white, skin, etc.
Packing: Simple packing, cutter, plastic can, Tinplate can



AL-WD30

Orthopedic Wool Padding

Material: 100% Cotton or Synthetic fiber
Weight: 390g/m², 400g/m², 410g/m², etc.
Size: 5cmx2.7m, 7.5cmx2.7m, 10cmx42.7, etc.
Packing: 1 roll / poly bag



AL-WD34

Non-woven Surgical Tape

Material: Non-woven with adhesive
Color: white, flesh
Packing: Simple packing, cutter, plastic can, Tinplate can



AL-WD35

PE Surgical Tape

Material: Polyethylene with adhesive
Color: transparent
Packing: Simple packing, cutter, plastic can, Tinplate can



AL-WD39

PU Adhesive Tape

Material: Polyurethane film
Color: Transparent
Size: 5cm×10m, 10cm×10m, 15cm×10cm, etc.
Packing: 100rolls/ctn, 50rolls/ctn



AL-WD36

Silk Surgical Tape

Material: silk with adhesive
Color: white or skin color
Packing: Simple packing, cutter, plastic can, Tinplate can



AL-WD40

Wound Adhesive Plaster

Cotton pad with/without medicine
Material: Cotton / PE / PU / PVC / Mesh / Cartoon / Spunlace, etc.
Size: Ø22mm, Ø25mm, 40×40mm, 38×38mm, 60×19mm, 72×19mm, etc.
Packing: 50pcs/box, 100pcs/box



AL-WD37

Fiberglass Casting Tape

Material: Fiberglass
Size: 5cm×360cm, 7.5cm×360cm, 10cm×360cm, etc.
Packing: 1pc/Aluminum foil bag, 50pcs/box



AL-WD41

Drilled Plaster

Material: Drilled zinc oxide with hot melt adhesive
Color: white, flesh
Size: 2.5cm×1m, 5cm×1m, 10cm×1m, etc.



AL-WD38

Paraffin Impregnated Gauze Pad

Sterile Paraffined Gauze Dressings are made of 100% cotton gauze swabs impregnated with white soft paraffin.
Size: 10cm×10cm, etc.



AL-WD42

Capsicum Plaster

Color: white, skin color
Size: 4cm×6cm, 7cm×10cm, 11cm×17cm, 12cm×18cm, etc.



AL-WD43

Self-adhesive with Absorbent Pad

Size: 6×3cm, 7×3.5cm, 7×5cm, 8×6cm, 10×10cm, 10×35cm, etc.

Packing: 1pc/blister pouch, 100pcs/box or 50pcs/box



AL-WD44

Waterproof and Transparent with Absorbent Pad

Size: 6×3cm, 7×3.5cm, 7×5cm, 8×6cm, 10×10cm, 10×35cm, etc.

Packing: 1pc/blister pouch, 100pcs/box or 50pcs/box



AL-WD45

Cannula Fixation Dressing (U shape)

Non-woven fabric material

Size: 7.6×5.1cm, 6×8cm, etc.

Packing: 1pc/blister pouch, 100pcs/box



AL-WD46

I.V. Cannula Dressing

Size: 6cm×7cm, 4.4cm×4.4cm, 7cm×8.5cm, etc.

Packing: 1pc/blister pouch, 100pcs/box



AL-WD47

Transparent Film Rolls

Material: PET backing, PU film, Medical pressure sensitive adhesive, Release paper

Size: 5cm×10m, 10cm×10m, 15cm×10m, etc.



AL-WD49

Adhesive Tape Remover Pad

Contents: Adhesive tape remover pad

Non-woven size: 65×30mm

Packing: 1pc/airtight pack, 100pcs/box



AL-WD51 Lemon Glycerine Swab Stick

Material: Impregnated with solution of water, Glycerine citric acid lemon flavouring, sodium, benzoate

Size: 70mm, 100mm, 150mm, etc.

Packing: 1pc or 2pcs or 3pcs/aluminum foil pouch



AL-WD48 Alcohol Prep Pad

Material: 70% Isopropyl Alcohol (75% Alcohol), Chlorhexidine

Inner Non-woven size: 65×30mm, 65×58mm, 80×40mm, 100×60mm, etc.

Packing: 1pc/airtight pack, 100pcs/box



AL-WD50

Cleaning Wipe (BZK)

Contents: Benzalkonium chloride or alcohol

Non-woven size: 200×130mm, 200×110mm, 45×45mm,

Packing: 110 wipes



AL-WD52

Alcohol Swab Stick

Material: 70% Isopropyl alcohol (75% alcohol)

Size: 70mm, 100mm, 150mm, etc.

Packing: 1pc or 2pcs or 3pcs/aluminum foil pouch



AL-WD53

Povidone-iodine Prep Pad

Material: 10% povidone-iodine solution, equivalent to 1% Available Iodine
Non-woven size: 65×30mm, packet size: 50×50mm
Packing: 1pc/airtight pack, 100pcs/box



AL-WD57

Dressing Set

1.1pc Sterile field
 2.5pcs Gauze swabs 7.5cm×7.5cm-12ply
 3.1pc Dressing tray
 4.1pc PE film 50×50cm
 5.1pc Yellow limpet bag 33×44cm
 6.1pc 10ml Syringe
 7.2pcs plastic forceps



AL-WD54

Povidone-iodine Prep Stick

Material: 10% povidone-iodine solution, equivalent to 1% Available Iodine
Size: 70mm, 100mm, 150mm, etc
Packing: 1pc or 2pcs or 3pcs/aluminum foil pouch



AL-WD58

Plastic Forceps

Material: PP, PE or ABS
Color: Blue, green, yellow, etc.
Size: 12cm, 14cm, 15cm, 19cm etc.



AL-WD55

Cleansing Wipe Tube

Contents: Benzalkonium chloride or Alcohol
 Inner Non-Woven
Size: 200×130mm, 110wipes
 200×110mm, 110wipes
 45×45mm, 110wipes



AL-WD59

Sponge Clean Stick

Size: 1#: 18cm flat stick, 2#: 15cm flat stick, 3#: 10cm circular stick
Packing: 20pcs/bag or 50pcs/bag



AL-WD56

Sting Relief

Inner Non-woven size: 65×30mm, 65×58mm, 80×40mm, 100×60mm, etc.
Packing: 1pc/airtight pack, 100pcs/box



AL-WD60

First-aid Box

1.1pc Quality PVC kit
 2.1bottle Polylab Antiseptic & Disinfectant Cleaning solution 60ml
 3.1bottle Polylab Povidone Iodine 10% w/v Broad Spectrum Antiseptic solution 15ml
 4.10pcs Fabric plaster strips
 5.10pcs Cotton balls
 6.100pcs Gauze swabs 5cm×5cm-8ply
 7.100tips Cotton buds
 8.1pc Elastic Crepe bandage 5cm×4.5m
 9.1pc W.O.W Bandage 7.5cm×4m
 10.1pc Triangular Bandage
 11.2pcs Gloves
 12.1roll Surgical tape
 13.1pc Scissors
 14.1bottle Medical oil 3ml
 15.1pc Usage card
 16.1pc First aid book



**Non-woven
& PPE**



AL-NW01

Round Nurse Cap

Material: Non-woven PP, 10g, 12g, 14g, 16g, etc.
Size: 18", 21", 24", etc.
Color: white, green, pink, blue, etc.
Packing: 100pcs/poly bag



AL-NW02

Clip Cap

Material: Non-woven PP, 10g, 12g, 14g, 16g, etc.
Size: 18", 21", 24", etc.
Type: Single elastic, double elastic
Color: white, green, pink, blue, etc.
Packing: 100pcs/poly bag



AL-NW03

Net Cap

Material: Nylon net fabric
Size: 18", 21", 24", etc.
Color: Blue, white, etc.
Packing: 100pcs/poly bag



AL-NW04

Doctor Cap

Material: Non-woven PP, 20g, 25g, 30g, 40g, etc.
Size: 61x14cm, etc.
Type: Tie back, Elastic back
Color: green, blue, etc.
Packing: 100pcs/poly bag



AL-NW06

Non-woven Face Mask

Ear-loop Type: 1ply, 2-plys or 3-plys
Ties on Type: 1ply, 2-plys or 3-plys
Color: blue, green, white and so on
Size: 17.5x9.5cm, etc.
Packing: 50pcs/box, 40boxes/ctn



AL-NW07

Nurse Paper Face Mask

2-ply or 1ply
Size: 20x7cm
Color: white
Packing: 100pcs/box, 100boxes/ctn



AL-NW05

Tudung/Hood Cap

Material: Non-woven PP, 20g, 25g, 30g, 40g, etc.
Size: 68×54cm, etc.
Type: Elastic
Color: green, blue, etc.
Packing: 20pcs/poly bag



AL-NW08

N95 Face Mask

With or without valve
Packing: 20pcs/box, 20boxes/ctn



AL-NW09

Beard Mask (Cover)

Specification:
Color: white, blue, green, etc.
Material: non-woven PP fabric
Gram: 10G / 12G / 14G / 16G / 18G / 20G etc.
Standard Size: 48×24cm



AL-NW10

Surgical Hood

Non Woven PP Material: 16g, 20g, 25g, 30g, etc.
Size: 30×124cm
Color: white, green, blue
Packing: 100pcs/poly bag, 10bags/carton



AL-NW12

Disposable Surgical Gown

Material: Non-woven PP, SMS, SMMS and so on
Color: white, blue, yellow, green, etc.
Weight: 40 / 45 / 50 / 55g etc.
Size: S, M, L, XL, XXL, XXXL
Packing: 1pc/sterilization blister bag, 50pcs/ctn



AL-NW13

CPE Gown

Material: CPE
Color: blue
Size: 125×200cm, etc.
Thickness: 0.03mm, etc.
Package: 1pc/bag, 100pcs/ctn



AL-NW14

Lab Coat

Material: non-woven PP/ PE
Color: white, blue, yellow, green and so on
Weight: 30/35/40/45g and so on
Size: M, L, XL, XXL
Packing: 10pcs/bag, 10bags/ctn



AL-NW15

Disposable Scrub Coat

Material: Non-woven PP, SMS
Weight: 48g, 55g and so on
Color: Blue, green and so on.
Packing: 1pc/bag, 50pcs/ctn



AL-NW16

Overall Gown

Material: PP, PE film, SMS or and so on with head hood, shoe cover and zipper
Color: white, blue, green, and so on
Weight: 35 / 45 / 50 / 60g and so on
Size: S, M, L, XL, XXL, XXXL
Packing: 1pc/bag, 50pcs/ctn



AL-NW11

Disposable Isolation Gown

Material: Non-woven PP+PE and so on
Color: white, blue, yellow, green and so on
Weight: 18 / 25 / 30 / 40g and so on
Size: S, M, L, XL, XXL, XXXL
Elastic cuff or knitting cuff
Packing: 10pcs/ bag, 100pcs/ctn





AL-NW17

4-ply Non-penetrated Anti-fog Face Mask

Specification:

Size: 9.2 (18)×18 cm, with anti-fog piece
Material: 18g/m² wet-laid nonwovens (or 22g/m² PP) + non-penetrated punching membrane+24g/m² melting-blown PP + 25g/m² calendering PP
Packing: 25pcs/box, 10box/carton
Function: 4-ply structure, effectively prevent from bacteria and dust, protect eyes and face from fluid, provide better safe property.
Production Technique: ultrasonic machine-sewed, with tie and antifog nosepiece
Application: suitable for casualty wards, infection wards, and scientific experiment



AL-NW18

Face Shield

Specification:

Size: 9.2 (18)×18 cm, with anti-fog piece
Material: 18g/m² wet-laid nonwovens (or 22g/m² PP) + non-penetrated punching membrane + 24g/m² melting-blown PP + 25g/m² calendering PP
Packing: 25pcs/box, 10box/carton
Function: 4-ply structure, effectively prevent from bacteria and dust, protect eyes and face from fluid, provide better safe property.
Production Technique: ultrasonic machine-sewed, with tie and antifog nosepiece
Application: suitable for casualty wards, infection wards, and scientific experiment



AL-NW19

Face Shield

Specification:

Efficiently two sides anti-fog film, high transparency. Excellent protection against potential contamination from blood pathogens, body fluids, or harmful chemical splash.
 Ventilation designed sponge, increase air flow and comfortable. Good bio-comparability, soft and comfortable to wear.
 Light and no pressure: wearing comfortable, with adjustable elastic headband.
 With excellent high and low temperature resistance.
Material: PET
Size: 32×22cm



AL-NW20

PVC Protective Goggles

Size: 15.0×7.0cm
Lens Material: PC/Optical PC
Frame Material: PVC/Medical PVC
Implementation Standard: GB14866-2006, EN166-2001, ANSI/ISEA Z87.1-2015
N.W./G.W.: 7.75KG/8.85KG
Carton Size A: 41.0×36.5×28.5CM, 60PCS/CTN
N.W./G.W.: 13.3KG/14.4KG
Carton Size B: 60.0×41.0×35.0CM, 200PCS/CTN



AL-NW22

PE Shoe Cover

Material: PE, 1.6g, 1.8g, 2.0g, 2.1g, 2.5g and so on
Size: 36×15cm, 41×15cm, 40×16cm and so on
Thickness: 0.016mm, 0.018mm, 0.02mm and so on
Packing: 100pcs/bag, 50bags/ctn



AL-NW21

Silicone Protective Goggles

Size: 15.0×7.0cm
Lens Material: PC/Optical PC
Frame Material: Medical grade silicone
Implementation Standard: GB14866-2006, EN166-2001, ANSI/ISEA Z87.1-2015
N.W./G.W.: 6.0KG/7.1KG
Carton Size: 41.0×36.5×28.5CM, 60PCS/CTN



AL-NW23


CPE Shoe Cover

Material: CPE, 1.8g, 2.0g, 2.1g, 2.5g and so on
Size: 36×15cm, 41×15cm, 40×16cm and so on
Thickness: 0.02mm, 0.025mm, 0.028mm and so on
Packing: 100pcs/bag, 50bags/ctn

AL-NW24

Non-woven Shoe Cover

Material: Non-woven
Color: blue, green, white and so on
Type A: PP regular shoe cover
Type B: PP Non-skid shoe cover
Weight: 30/40/45g and so on
Size: 36×15cm, 38×18cm, 40×16cm and so on
Packing: 100pcs/bag, 10bags/ctn





AL-NW25

Boot Cover


Material: Breathable Non-woven 63g, etc.
Color: White color with blue belt
Size: 40cm×43cm or another can be customized
Packing: 50pcs/bag, 500pcs/ctn



AL-NW26

Non-woven Boot Cover

Material: Non-woven material 40g, etc.
Color: Blue, green, etc.
Size: 36×54cm or another can be customized
Package: 50pcs/bag, 500pcs/ctn



AL-NW31

Non-woven Bed Cover

Size: 210×90×20cm and so on
Weight: 30g and so on
Color: blue, green, white and so on
Packing: 1pc/roll, 10rolls/bag, 2bags/ctn



AL-NW32

Non-woven Pillow Cover

Size: 50×60cm with 5cm gusset, and so on
Weight: 80g and so on
Color: blue, green, white and so on
Packing: 100pcs/bag, 1bag/ctn



AL-NW27

PE Boot Cover

Material: 10g/pc, etc.
Color: Blue, white, etc.
Size: 36×45cm or another can be customized
Package: 50pcs/bag, 500pcs/ctn



AL-NW28

PE Sleeve

Size: 20×40cm, 42×22cm, 46×22cm and so on.
Thickness: 0.02mm, 0.018mm
Color: white, blue, red and so on
Packing: 100pcs/bag, 20bags/ctn



AL-NW33

Diaper

Type: Adult / Child
Material: Inner layer: Non-woven
Absorption: Cotton pulp, Tissue, SAP
Stopper: Adhesive tape, Frontal tape
Back Sheet: water-proof PE back layer
Other: Leak guard, Fabric glue and so on
Size: S / M / L / XL
Packing: 10pcs/bag, 8 bags/ctn



AL-NW34

Disposable Underpad

Material: Non-woven fabrics+pulp+sap+PE film
Color: white, blue, green and so on
Size: 60×60cm, 60×90cm and so on
Packing: 10pcs/bag, 10bags/ctn



AL-NW29

PE Apron

Material: PE, with adjustable Tie
Size: 110×70cm, 135×71cm 160×90cm and so on
Thickness: 0.018mm, 0.020mm, 0.025mm
Color: white, blue and so on
Package: 10pcs/bag, 100bags/ctn



AL-NW30

Non-woven Bed Sheet

Size: 60×180cm and so on
Weight: 30g and so on
Color: blue, green, white and so on
Packing: 25pcs/poly bag, 40bags/ctn



AL-NW35

Body Bag

With handle or without handle
Color: White, Black, etc.
Material: PEVA, PVC, Nylon, etc.
Size: Adult 230×92cm, Child 120×80cm, etc.

Home Care



AL-HC01

PVC Hot Water Bag

Material: PVC
Size: 750ml-1000ml-1750ml-2000ml-2500ml-3000ml
Color: blue, green, red



AL-HC02

Latex Hot Water Bag

Material: Latex
Size: 750ml-1000ml-1750ml-2000ml-2500ml-3000ml
Color: blue, green, red



AL-HC03

Cover for Hot Water Bag

Material: flannelette/knit/furry etc
Size: 750ml-1000ml-1750ml-2000ml-2500ml-3000ml



AL-HC04

Ice Bag

Size: 6"/9"/11"



AL-HC05

Instant Ice Pack

Size: 280g, 23×14CM
160g, 19.5×12CM
100g, 15×11CM



AL-HC06

Instant Hot & Cold

Size: 250g 23×13cm, 70g 11×9cm, 70g 11×10.5cm



AL-HC07 Cold Eye Mask

Size: 30g 20×5cm, 50g 23×7cm, 140g 41×10cm



AL-HC08 Air Cushion Latex

Material: calcium carbonate+rubber
Size: 35cm / 40cm / 42.5cm / 45cm



AL-HC09 Ear Syringe

PVC, 30ML / 60ML / 90ML



AL-HC10 Ear Syringe

Latex, 30ML / 60ML / 90ML



AL-HC11

Cooling Gel Patch

For instant cooling sensation soothes fever discomfort
Produce fever quickly-easy to use
Gentle to the skin
Safe to use with medications
Contain non-drug
Soothes fever instantly
Strong adhesive power
10 hours cooling effect

Use Instructions:

Pull off the transparent film and apply aikido gel surface to the effeled area.
to enhance the cooling effect, place aikido gel in a refrigerator before use, don't place aikido gel in freezer. Aikido gel may be cut to the required size using a pair of scissors



AL-HC12

Thermometer Cover

Material: PE
Description: Insert, twist, pull apart or peel back covers
Packing: 200pcs/box, 150boxes/carton
Carton Dimension: 38×32×32cm



AL-HC13

I.D.Bracelet

Usage: Record Important Identification Information of A patient, Tie on patient's wrist, 24-hour close identification can effectively ensure that any time quickly and accurately identify patients.

Material: PVC
Color: Blue, Pink, White, Other color by optinal
Model: Adhesive, write-on, Insert-card, infant, Adult, Adult & Infant

AL-HC00

Instant Hand Sanitizer

Volume: 0.5ml/15ml/29ml/53ml/59ml/236ml/500ml/1000ml




AL-HC14

Ear Thermometer Probe Cover



AL-HC15 4 Case Pill Box



AL-HC17 7 Case Pill Box



AL-HC19 14 Case Pill Box



AL-HC21 Pill Crusher

Material: ABS+PS, Feature Crush and Store Pills
Colour: Blue, Red, White, Function easy to use for pill Storage and crushing



AL-HC16 6 Case Pill Box



AL-HC18 8 Case Pill Box




AL-HC20 28 Case Pill Box



AL-HC22 Pill Cutter

Material: PS
Size: 9.5×4.5×3cm



AL-HC23

Heat-Sealing Gusseted (Reel)

Size:
 (75+25)mm×100mm, 1roll/box, 5rolls/ctn,
 22.5×46×23cm
 (100+50)mm×100mm, 1roll/box, 4rolls/ctn,
 22.5×47×23cm
 (150+50)mm×100mm, 1roll/box, 2rolls/ctn,
 22.5×34×23cm
 (400+75)mm×100mm, non-box, 1rolls/ctn,
 22×41×22cm



AL-HC27

Heat-Sealing Flat Reel Pouch

Size:
 55mm×200mm, 1roll/box, 6rolls/ctn, 22.5×42×23cm
 75mm×200mm, 1roll/box, 5rolls/ctn, 22.5×46×23cm
 100mm×200mm, 1roll/box, 4rolls/ctn, 22.5×47×23cm
 150mm×200mm, 1roll/box, 2rolls/ctn, 22.5×34×23cm
 400mm×200mm, non-box, 1rolls/ctn, 22×41.5×22cm



AL-HC24

Self-sealing Flat Pouch

Size:
 50×130mm 200pcs/box, 14400pcs/ctn,
 270.6×37×30cm
 90×230mm 200pcs/box, 4800pcs/ctn,
 38.2×23×30cm
 135×255mm 200pcs/box, 3600pcs/ctn,
 26.8×42.4×30cm
 190×330mm 200pcs/box, 2400pcs/ctn,
 34.3×29.6×30cm
 250×370mm 200pcs/box, 1200pcs/ctn,
 38.3×26.3×30cm



AL-HC28

Medical Crepe Paper

Size:
 60×60mm
 75×75mm
 100×100mm, etc.



AL-HC25

Paper/Paper Gusseted Pouch


Size:
 90×50×150mm
 90×50×250mm
 110×30×190mm
 140×75×250mm, etc.



AL-HC29

Sterilization Indicator Tape

Size: 12.5mm×50m, 19mm×50m, 25mm×50m etc.
Color: Light Yellow
Color range: Light yellow→Black



AL-HC26

Heat-Sealing Gusseted Pouch

Size:
 (75+25)mm×300mm, 100pcs/box, 3000pcs/box,
 31.3×40×30cm
 (100+50)mm×300mm, 100pcs/box, 2400pcs/box,
 31.3×42.2×30cm
 (150+50)mm×400mm, 100pcs/box, 1200pcs/box,
 41.3×31.6×30cm
 (200+55)mm×400mm, 100pcs/box, 600pcs/box,
 21.3×41.3×30cm
 (250+60)mm×450mm, 100pcs/box, 600pcs/box,
 26.3×46.3×30cm



AL-HC30

Steam Indicator Strip


Size: 15×100mm, 15×200mm, 200pcs/box, etc.
Color range: Blue→Black



ECG

ECG Paper



Single-lead ECG Paper
3-lead ECG Paper
6-lead ECG Paper
12-lead ECG Paper



210mmX300mm-200P


12-lead ECG Paper

AL-EC01
ECG Paper




50mmX20m-30m

Single-lead ECG Paper



6-lead ECG Paper



3-lead ECG Paper

Single-lead ECG Paper: 50mm×20m, 200pcs/ctn
3-lead ECG Paper: 90mm×90mm-200P, 100pcs/ctn
6-lead ECG Paper: 112mm×90mm-200P, 100pcs/ctn
12-lead ECG Paper: 210mm×300mm-200P, 20pcs/ctn



AL-EC02

Disposable Oxi-Sensor

9 pin, Neonate/Adult
1pc/blister packing, 40pcs/box, 4boxes/ctn



Disposable Electrode

To comply with CE, ISO 13485:2003 and FDA.
Comfortable and Convenient, no residual and non-irritable
Excellent adhesive from special gel formula
Different sizes and material have been studied to meet the different requirements.

General Monitoring

Foam Recommended for ICU/CCU, Surgery, Emergency

Item Code	Description	Packaging	Case Dimension
AL-EC03-1	Child, Solid Gel, 36x44mm	50/pk, 40pk/cs	48x36x28cm
AL-EC03-2	Adult, Solid Gel, 42x50mm	50/pk, 40pk/cs	48x36x28cm
AL-EC03-3	Adult, Solid Gel, Ø50mm	50/pk, 40pk/cs	48x36x28cm
AL-EC03-4	With 4mm banana plug, Ø50mm	50/pk, 40pk/cs	43x55x28cm
AL-EC03-5	Neonatal, Solid Gel, Ø25.4mm	50/pk, 40pk/cs	40x36x28cm

Cloth Recommended for ICU/CCU, Emergency, Long term

Item Code	Description	Packaging	Case Dimension
AL-EC04-1	Cloth, Solid Gel, Ø50mm	50/pk, 40pk/cs	40x36x28cm

Clear Tape Recommended for ICU/CCU, Surgery, Emergency, Long term, Holter

Item Code	Description	Packaging	Case Dimension
AL-EC05-1	Breathable clear tape, Solid Gel, Ø50mm	50/pk, 40pk/cs	40x36x28cm



AL-EC06

Medical Conductive Paste

Property unit value
sound speed(35°C)M/s 1520-1620
Electric resistance(35°C)-13-17
Unit area acoustic impedance (35°C)
Pa-s/m 1.5-1.7×16 sound attenuation
(35°C)Db/(cm.mhz)≤0.05
Viscosity (25°C) Pa.s ≥15
Ph-5.5-8.0
Size: 60ml, 100ml, 150ml, 250ml



AL-EC07



AL-EC08

AL-EC07/08

Ultrasonic Gel

Property unit value
sound speed(35°C)M/s 1520-1620
Unit area acoustic impedance (35°C)
Pa-s/m 1.5-1.7×16 sound attenuation
(35°C)Db/(cm.mhz)≤0.05
Viscosity (25°C) Pa.s ≥15
Ph-5.5-8.0
Size: 250ml, 1000ml, 2500ml, 5000ml



Radilo Translucent

Foam Recommended for Surgery, Emergency, Cath Lab, Radiology

Item Code	Description	Packaging	Case Dimension
AL-EC08-1	Adult, Solid Gel, Ø50mm	50/pk, 40pk/cs	40x36x28cm

Stress / Diaphoretic

Recommended for ICU/CCU, Surgery, Stress, Holter, Emergency

Item Code	Description	Packaging	Case Dimension
AL-EC09-1	Adult, Solid Gel, 42x37mm	50/pk, 40pk/cs	40x36x28cm

Wet Gel


Recommended usage ICU/CCU, Surgery, Stress, Holter, EMS, Emergency

Item Code	Description	Packaging	Case Dimension
AL-EC10-1	Adult, Wet Gel, 52x47mm	50/pk, 40pk/cs	40x36x28cm

Repositionable

Recommended usage ICU/CCU, Surgery, Emergency, Cath Lab, Radiology

Item Code	Description	Packaging	Case Dimension
AL-EC11-1	Adult, Solid Gel, 32x40mm	3/pk, 10pk/bx, 12bx/cs	48x22x22cm



AL-EC12

ECG Resting Electrodes

Size:
14×30mm
20×30mm
26×30mm
23×34mm
28×34mm

Diagnostics

Stethoscope

1. **Material of chest piece:**
Zinc Alloy, aluminum Alloy, Stainless steel
2. **Material of earloop:** Copper, Aluminum, Stainless steel.
3. **Material of earplug:** PVC, ABS
4. **Material of tube:** PVC.
5. Any color of tube is available.



AL-DG01

Pediatric
single head stethoscope



AL-DG02

Deluxe Special

Single head stethoscope
Zinc alloy chest piece
with chrome plated silver finish



AL-DG03

Dual head stethoscope
Colored chest piece & Binaural

Sprague Rappaport Stethoscope

Sprague Rappaport type chest piece
Brass chrome plated binaural
Plastic or soft eartips
All kinds of color available for the such as:
Black, Red, Green, Yellow, Blue, Purple,
Pink, Grey, White, etc.
2sets of earplug, 1set of mebrane
1pc of large chest piece
1pc of middle chest piece
1pc of small chest piece



AL-DG04

Sprague Rappaport Stethoscope
Colored chest piece

Stainless steel stethoscope

Stainless steel chest piece
Stainless steel Binaural
Plastic & Soft earrips
With non-chill rings
Spare diaphragm included
All kinds of color available for the such as:
Black, Red, Green, Yellow, Blue, Purple,
Pink, Grey, White, etc.



AL-DG05

Stainless steel stethoscope



AL-DG09

Aneroid Sphygmomanometer



AL-DG10

Desk type
Aneroid Sphygmomanometer



AL-DG06

Stainless steel stethoscope
Adult type single-lumen tubing



AL-DG07

Stainless steel stethoscope
Child type



AL-DG11

Aneroid Sphygmomanometer with
separate single head stethoscope



AL-DG08

Master Cardiology Stethoscope

Specification:
Application: Cardiology/excellent performance.
Product weight/length: 215g/71cm
Chest piece material: stainless steel,
Surface treatment: manual polish
Diaphragm: $\phi 44\text{mm}$, **Type:** floating diaphragm
Binaural material: stainless steel, $\phi 6\text{mm}$
Eartip material: silicone
Structure: Soft sealing structure
Accessory box: 2 sets silicone ear tips
Package Content: 1pc/box, 12pcs/ctn,
580x400x245mm



AL-DG12

Wall type
Aneroid Sphygmomanometer



AL-DG13

Floor Type Aneroid Sphygmomanometer

Form of packing:
One machine fill in a white box, and every one boxes
in one export carton.
Size of packing:
White box: 495x460x170mm
Export carton: 500x470x185mm



Mercury Sphygmomanometer

300mmHg scale glass tubes Metal case
Cotton or Nylon cuff
With or without D-ring is available
Artery Label
Latex Bladder with double tubes



AL-DG14
LCD Display



AL-DG15
Light Display

Accessories



AL-DG17
Zipper Case



AL-DG16
Bulb



AL-DG18
Bladder



AL-DG20
Cuff



AL-DG19
Air-release Valve



AL-DG21
End Valve



AL-DG22

Blood Pressure Monitor

Measurement method: oscillographic method
Accuracy: $\pm 3\text{mmHg} \pm 0.4\text{kPa}$ (blood pressure), $\pm 5\%$ (pulse)
Measurement range: 0-290mmHg (0-39kPa)
Power supply: lithium battery DC3.7V
Dimension: 137×109×69mm
Weight: 250g with battery
Memory recall: 2×90 groups
Auto power-off : 60s no-operation



AL-DG23

Blood Pressure Monitor

Measurement method: oscillographic method
Accuracy: $\pm 3\text{mmHg} \pm 0.4\text{kPa}$ (blood pressure), $\pm 5\%$ (pulse)
Measurement range: 0-290mmHg (0-39kPa)
Power supply: 4XAAA batteries (DC6V)
Dimension: 137mm×109mm×71mm
Weight: 256g with battery
Memory recall: 2×90 groups
Auto power-off : 60s no-operation



AL-DG26

Fingertip Pulse Oximeter

Features:
 Parameter: SPO₂, HR, PI
 TFT color screen with four direction display
 One button operation, simple and easy
 8s auto power off without operation
 Compact body and portable to carry
Specification:
 Body weight: 44g (Without batteries)
 Dimensions: 62×31×36mm
 Security type: class 2 BF type application equipment
 Screen type: LED/TFT/OLED color display
 Product Specification: finger clip type
 Battery: 2×AAA 1.5V alkaline batteries
 SpO₂ measuring range: 0-100% (the resolution is 1%)
 PR measuring range: 30-250BPM (the resolution is 1BPM)
 SpO₂ measuring accuracy: 70-100%; $\pm 2\%$, $\leq 70\%$ unspecified
 PR measuring accuracy: $\pm 1\text{BPM}$ or $\pm 1\%$ (select larger)
 Package: Oximeter, user manual, warranty card



AL-DG27

Fingertip Pulse Oximeter

Features:
 Parameter: SPO₂, HR
 LED color screen with four direction display
 One button operation, simple and easy
 8s auto power off without operation
 Compact body and portable to carry
Specification:
 Body weight: 43g (Without batteries)
 Dimensions: 62×31×36mm
 Security type: class 2 BF type application equipment
 Screen type: LED/TFT/OLED color display
 Product Specification: finger clip type
 Battery: 2×AAA 1.5V alkaline batteries
 SpO₂ measuring range: 0-100% (the resolution is 1%)
 PR measuring range: 30-250BPM (the resolution is 1BPM)
 SpO₂ measuring accuracy: 70-100%; $\pm 2\%$, $\leq 70\%$ unspecified
 PR measuring accuracy: $\pm 1\text{BPM}$ or $\pm 1\%$ (select larger)
 Package: Oximeter, user manual, warranty card



AL-DG24

Blood Pressure Monitor

Measurement method: oscillographic method
Accuracy: $\pm 3\text{mmHg} \pm 0.4\text{kPa}$ (blood pressure), $\pm 5\%$ (pulse)
Measurement range: 0-290mmHg (0-39kPa)
Power supply: lithium battery DC3.7V
Dimension: 137×109×69mm
Weight: 250g with battery
Memory recall: 2×90 groups
Auto power-off : 60s no-operation



AL-DG25

Blood Pressure Monitor

Measurement method: oscillographic method
Accuracy: $\pm 3\text{mmHg} \pm 0.4\text{kPa}$ (blood pressure), $\pm 5\%$ (pulse)
Measurement range: 0-290mmHg (0-39kPa)
Power supply: 4XAAA batteries (DC6V)
Dimension: 114mm×89mm×126mm
Weight: 264g with battery
Memory recall: 2×90 groups
Auto power-off : 60s no-operation



AL-DG28

Fingertip Pulse Oximeter

Specification:
 Body weight: 41.5g (Without batteries)
 Dimensions: 59×31×29mm
 Security type: class 2 BF type application equipment
 Screen type: TFT/LED color display
 Product Specification: finger clip type
 Battery: 2×AAA 1.5V alkaline batteries
 SpO₂ measuring range: 0-100% (the resolution is 1%)
 PR measuring range: 30-250BPM (the resolution is 1BPM)
 SpO₂ measuring accuracy: 70-100%; $\pm 2\%$, $\leq 70\%$ unspecified
 PR measuring accuracy: $\pm 1\text{BPM}$ or $\pm 1\%$ (select larger)
 Package: Oximeter, user manual, warranty card



AL-DG29

Fingertip Pulse Oximeter

Specification:
 Body weight: 41.5g (Without batteries)
 Dimensions: 58×31×32mm
 Security type: class 2 BF type application equipment
 Screen type: TFT/LED color display
 Product Specification: finger clip type
 Battery: 2×AAA 1.5V alkaline batteries
 SpO₂ measuring range: 0-100% (the resolution is 1%)
 PR measuring range: 30-250BPM (the resolution is 1BPM)
 SpO₂ measuring accuracy: 70-100%; $\pm 2\%$, $\leq 70\%$ unspecified
 PR measuring accuracy: $\pm 1\text{BPM}$ or $\pm 1\%$ (select larger)
 Package: Oximeter, user manual, warranty card



AL-DG30

Digital Thermometer

Large LCD Display
Type: Flexible Tip
Display range: 32.0°C - 42.9°C (90.0°F - 109.9°F)
Accuracy: $\pm 0.1^\circ\text{C}$ ($\pm 0.2^\circ\text{F}$)
Min. scale: 0.1
Measurement time: 60s
Beeper, Auto Shut-off, Memory
Optional: waterproof



AL-DG31

Digital Thermometer

Large LCD Display
Type: Standard Tip
Display range: 32.0°C - 42.9°C (90.0°F - 109.9°F)
Accuracy: $\pm 0.1^\circ\text{C}$ ($\pm 0.2^\circ\text{F}$)
Min. scale: 0.1
Measurement time: 60s
Beeper, Auto Shut-off, Memory
Optional: waterproof



AL-DG32

Infrared Thermometer

Measuring range: 32.0°C~42.2°C
Measuring mode: forehead mode
Memory recall: 32 memories recall
Resolution: 0.1°C
Dimension: 146mm×89mm×40mm
Power supply: 2×AAA batteries
(DC 3V), replaceable
Weight: 138g without battery



AL-DG33

Infrared Thermometer

Measuring range: 32.0°C~42.2°C
Measuring mode: forehead mode
Memory recall: 32 memories recall
Resolution: 0.1°C
Dimension: 144mm×101mm×26mm
Power supply: 2×AAA batteries
(DC 3V), replaceable
Weight: 88g without battery



AL-DG34

Infrared Thermometer

Measuring range: 32.0°C~42.2°C
Measuring mode: forehead mode
Memory recall: 32 memories recall
Resolution: 0.1°C
Dimension: 141mm×47mm×37mm
Power supply: 2×AAA batteries
(DC 3V), C/F conversion
Weight: 64g without battery



AL-DG35

Infrared Thermometer

Measuring range: 32.0°C~42.2°C
Measuring mode: forehead mode
Memory recall: 32 memories recall
Resolution: 0.1°C
Dimension: 138mm×47mm×37mm
Power supply: 2×AAA batteries
(DC 3V), replaceable
Weight: 47g without battery



AL-DG36

Infrared Thermometer

Measuring range: 32.0°C~42.2°C
Measuring mode: forehead mode
Memory recall: 32 memories recall
Resolution: 0.1°C
Dimension: 164mm×44mm×40mm
Power supply: 2×AAA batteries
(DC 3V), replaceable
Weight: 31g without battery



AL-DG37

Infrared Thermometer

3-Digital, Scale: °C or °F
Accuracy: $\pm 0.2^\circ\text{C}$ ($\pm 0.4^\circ\text{F}$)
Min Scale: 0.1 °C
Memory Function: 10 group
Measuring Times: 1 second
Beeper function, Auto shut-off



AL-DG38



AL-DG39



AL-DG40

AL-DG38

Buck Neurological Hammer

Hammerhead: zinc alloy with rubber plug
Handle: Nylon
With needle: Zinc alloy, with brush
Length of hammer: 22.5cm
Height of hammer head: 8.5cm

AL-DG39

Gear Wheel Hammer

Hammerhead: Zinc Alloy
Handle: Zinc Alloy

AL-DG40

Taylor Hammer

Length of Handle: 18cm
Size of Runner Plug: 5×4.3×1cm
Thickness of Rubber Plug: 1.5cm



AL-DG41

Gristletly Hammer

Rubber Plug
Hammerhead: copper
Handle: steel
Length of hammer: 24cm
Diameter of hammerhead: 4.5cm



AL-DG42

Babinski (Adult)

with Telescoping Handle
Hammerhead: Zinc alloy with rubber plug
Handle: Zinc Alloy
Length of handle: 21cm
Diameter of hammerhead: 4.5cm
Thickness: 1cm



AL-DG43

Colorflex Reflex Hammer

Hammerhead: rubber with copper ring
Handle: Nylon
Length of hammer: 22.5cm
Height of hammerhead: 7.5cm



AL-DG44

Diagnostic Pen Light

AL-DG45

Diagnostic Pen Light



AL-DG46 Otoscope Set

Specification:

1. One chargeable Otoscope & Ophthalmoscope
2. Single Charger
3. Gift box

AL-DG47 Otoscope Set

Specification:

1. One chargeable Otoscope & Ophthalmoscope
2. Double Charger
3. Gift box



AL-DG49

ENT Kit

Ophthalmoscope	1
Otoscope carrier	1
Ear specula (3, 4, 5, 6mm diam)	4
Nasal speculum (8mm diam)	
Dilating Nasal speculum	1
Laryngeal mirror	3
Tongue depressor plastic	1
Battery handle	1
Spare Bulbs (2large, 4small)	6



AL-DG50

Professional Otoscope

Features:

Fiberoptics for cool, obstruction free light transmission
Otoscope head made from polymid - an impact resistant space age thermoplastic
Swivelling viewing window with 2.5x magnification "AA" polycarbonate battery handle interchanges pocket clip with integral on/off slide switch for ambidextrous use
Technical Data

Magnification: 2.5x, **Distal dia.:** 2.5mm, 4mm

Bulb: 2.5V 1.6W halogen lamp

Batteries: two "AA" Batteries or two "AA" NiMH Rechargeable Battery



AL-DG51

Otoscope

Deluxe packing:
Penlight, Part of otoscope, batteries



AL-DG52

Reflex Hammer Gifts Set

Buck Neurological Hammer
Babinski with telescoping handle
Taylor Hammer



AL-DG53

Buck Neurological Hammer

Gear wheel hammer
Taylor Hammer
Gristlet hammer



AL-DG54

Fiber Laryngeal Mirror

Power supply: 1.5V X 2 "C" Type
Bulb: 2.5V / 0.35A
Blades: 92, 114, 134mm



AL-DG48

Integrated Diagnostic System

Fiber Optic Otoscope
Direct Ophthalmoscope
Tip Dispenser (with disposable tips)
Infrared thermometer
Insufflation Bulb
Sphygmomanometer
Customize your wall set to match your individual needs



AL-DG55

Height Measuring Rod

Height range to be measured: 0-2000mm
Min value of height per division: 1mm
Measurement: 35,5×19×19cm / 20pc
G.W.: 4kg, N.W.: 3.6kg



AL-DG57

Baby Scale

Max. Weighing: 20kg



AL-DG59

Hanging Scale with Panty and Bag

Max. Weighing: 25kg



AL-DG56

Electronic Infant Scale

Max. Weighing: 20kg



AL-DG58

Electronic Infant Scale

Max. Weight: 20kg
Min. Value per division: 50g



AL-DG60

Health Scale

Max. Weighing: 220kg



AL-DG61

Health Scale

Max. Weighing: 200kg

Medical Oxygen Regulator

Input Pressure: 15Mpa
output Pressure: 0.2-0.3Mpa
Automatic Discharge Pressure of The Safety Valve: 0.35±0.05Mpa
The Range of Flow: 1-15L/Min
Cylindercconnecting nut: G5/8 right hand (according to the need can be customized)
Piton decompression structure
Flow adjusting knob: adjust flow easily, shut off the flow without oxygen output
Tide change glass(bottle): 150ml, immerses disinfection, can be 121 °C high tempreature disinfection, bottles have water level indicator

Medium	Metering range (L/min)	Accuracy	Entrance Pressure (MPa)	Exit Pressure (MPa)	Entrance connection screw thread	Output connection screw thread
O ₂	1-15	grade 4	15MPa 20MPa	0.2-0.3 MPa	G5/8-14-RH CGA540-RH G3/4-14-RH W21.8-14-RH Pin-index	8 mm



AL-DG62



AL-DG63



AL-DG64



AL-DG65

Seamless Steel Cylinder, Aluminum cylinder (portable oxygen cylinder bag, oxygen cylinder case)



CAT NO.	Height (MM)	Weight (KG)	Water Capacity (L)	Nominal outer Dia. (MM)	Nominal working pressure (bar)	Water test pressure (bar)	Design Wall Thickness (MM)	Material	Production Standard	Remark
AL-89-1-15	237	1.76	1	89	150	250	3	37Mn	GB5099 ISO9809	Bottle
AL-108-2-15	310	3.6	2	108	150	250	4	37Mn	GB5099 ISO9809	Bottle
AL-108-3-15	437	5	3	108	150	250	4	37Mn	GB5099 ISO9809	Bottle
AL-108-4-15	565	6.2	4	108	150	250	4	37Mn	GB5099 ISO9809	Bottle
AL-140-4-15	365	6.6	4	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-5-15	440	7.7	5	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-6-15	515	8.8	6	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-7-15	595	9.9	7	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-8-15	665	11	8	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-9-15	745	12.2	9	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-10-15	830	13.4	10	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-140-13.4-15	1070	16.8	13.4	140	150	250	4.5	37Mn	GB5099 S09809	Bottle
AL-140-14-15	1115	17.7	14	140	150	250	4.5	37Mn	GB5099 ISO9809	Bottle
AL-159-10-15	665	14.1	10	159	150	250	5	37Mn	GB5099 1509809	Bottle
AL-159-12-15	790	16.5	12	159	150	250	5	37Mn	GB5099 ISO9809	Bottle
AL-159-14-15	890	18.3	14	159	150	250	5	37Mn	GB5099 S09809	Bottle
AL-159-15-15	945	20.1	15	159	150	250	5	37Mn	GB5099 ISO9809	Bottle
AL-18020-15	1010	29.2	20	180	150	250	6	37Mn	GB5099 ISO9809	Bottle
AL-219 40-15	1350	48.5	40	219	150	250	5.7	37Mn	GB5099 ISO9809	Bottle
AL-232-45-15	1327	53.8	45	232	150	250	6	37Mn	GB5099 ISO9809	Bottle
AL-232 47-15	1379	55.7	47	232	150	250	6	37Mn	GB5099 ISO9809	Bottle
AL-232-50-15	1457	58.5	50	232	150	250	6	37Mn	GB5099 ISO9809	Bottle



AL-DG66

Oxygen Concentrator

Product Size: 47×28.5×55.6cm

Oxygen flow: 0-3L/min, 0-5L/min, 0-8L/min, 0-10L/min

Power: 330W, 500W, 650W, 700W

Noise: 48dB(A), 52dB(A), 55dB(A), 59dB(A)

Net weight: 18.5kgs, 23kgs, 25kgs, 27kgs

Oxygen density: 93%±3% With nebulizer function

Nebulizer rate: 0.1mL/min

Optional function: oxygen density alarm function



AL-DG68

Infant Incubator

Specification

Power requirement: AC220V-230V/50Hz

or AC110-120V/50-60Hz, 850VA

Control Mode: Air mode control

Air temperature control range: 25℃-37℃(override mode: 37℃-39℃ is optional)

The temperature fluctuation range between incubator temperature and average incubator temperature (Under stable temperature condition): ≤0.5℃

Temperature uniformity: ≤0.8℃(mattress is on the horizontal position)

Inclination of the bassinet: ± 12°

Noise inside hood: ≤45dB(A)

Failure alarm: over temperature alarm, deviation alarm, sensor failure alarm, fan motor alarm, power failure alarm, system alarm and so on.

Built-in timer range: 0~9999 hrs 59min

Accuracy of weight display: ± 1% (equipped with weighing system)



AL-DG67

Oxygen Concentrator

Product Size: 43×31×56cm

Oxygen flow: 0-3L/min, 0-5L/min

Power: 330W, 500W

Noise: 48dB(A), 52dB(A)

Net weight: 18.5kgs, 21kgs

Oxygen density: 93%±3% With nebulizer function

Nebulizer rate: 0.1mL/min

Optional function: oxygen density alarm function



AL-DG69

Full Digital Ultrasound Scanner

Technical Parameters:

Image Storage: 64 Frames

Display Depth: ≥160

Scanning Model: Electronic Linear Array, Electronic Convex Array

LCD Size: 9.7 Inch/ 12Inch

Gray Scale: 256 Levels

Image Flip: Up/ Down, Left/ Right, Blake/ White

Image Processing: Image Smoothing/Sharpening, Tissue Harmonic, Histogram, DR, Gamma Correction, Pseudo-Color

Measurement: Freehand, Ellips, Distance for Perimeter, Area Volume, Heart, G.A, EDD, BPD-FW, FL, AC, HC, CRL, AD, GS, LMP

Character Display: Date, Clock, Name, PID, Age, Sex, Hospital Name, Doctor

Notation: Full-screen Character Editor, Body Mark, Position Indication

Report: Abdomen and gynecology, Cardiac, Obstetrics, Urology

Display Mode: B, B+B, B+M, M, 4B

Scanning Angle: Adjustable

Cine Loop: ≥400Frames

Output: VGA and PAL video outputs

Power Consumption MAX: 100VA

Product Size: 289×304×222mm

Carton Size: 395×300×410mm

N.W./ G.W.: 6kg/ 7kg

Tell us what you need.

Send us the product codes, specifications, expected volume, destination market and delivery timing. Our team typically responds within one business day.

Hong Kong

HEADQUARTERS

1505–8 Hing Yip Commercial Building
284 Des Voeux Road Central
Hong Kong

+852 9811 8857

Dubai

REGIONAL OFFICE

297 Bayan Business Center
Green Community Village
Dubai, UAE

+971 50 135 5145

EMAIL info@alitmedical.com

WEBSITE www.alitmedical.com



We source it.
We inspect it.
We deliver it.

www.alitmedical.com